

## Protocol

# Immunocytochemistry (ICC)

**Note:** This procedure is optional if detecting a membrane protein.

## Reagents Required

Product	Preparation	Item No.
Phosphate Buffered Saline (PBS)		MB-008
Triton X-100		
Hydrogen Peroxide (H <sub>2</sub> O <sub>2</sub> )		
Bovine Serum Albumin (BSA)		BSA-10
Biotinylated Secondary Antibody	Dilute 1:500	
Streptavidin Peroxidase Conjugated	Dilute 1:500	S000-03
3,3'-Diaminobenzidine (DAB) Substrate		DAB-10
Hematoxylin (Optional) and Acetic Acid (Optional) for counterstain		KHG001
UltraPure Sterile Water		MB-009-1000
Coverslip Solution	50% glycerol/UltraPure water	

## Procedure

1. Add one drop of PBS/0.1% Triton X-100 to each well to permeabilize the cells and incubate slides for one 1 minute at room temperature.
2. Remove the liquid and wash the slides twice in PBS, 5 minutes each on the shaker.
3. Remove the liquid and place the slides onto a tray.
4. Soak slides in 1.5% H<sub>2</sub>O<sub>2</sub>/PBS solution for 15 minutes.
5. Wash twice in PBS for 5 minutes each on the shaker.
6. Incubate with 5% BSA into each well to block overnight at 4°C in a humid chamber.
7. Dilute the primary antibody to the recommended concentration in 1% BSA diluent.
8. Remove BSA from the slides.
9. Add 35 µL of primary antibody to each well. Incubate for one 1 hour at room temperature.
10. Remove the primary antibody solution and wash slides 3 times in PBS, 5 minutes each on the shaker.
11. Dilute the biotinylated secondary antibody to 1:200 in a solution of 1% BSA diluent.

12. Remove the excess fluid and add one drop secondary antibody solution into each well. Incubate for one 1 hour at room temperature.
13. Wash in PBS 3 times, 5 minutes each on an orbital shaker. Remove excess fluid.
14. Add one drop streptavidin peroxidase to each well. Incubate for 30 minutes at room temperature.
15. Wash 3 times, 5 minutes in PBS on an orbital shaker. Remove excess fluid.
16. Add DAB substrate to each cell well. Once the cells start turning brown wash 2 times in PBS for 5 minutes each on the shaker.  
**Note:** Inexperienced technicians may wish to observe cells turning brown under a microscope.
17. (Optional step for counterstain) Dip the slide rack with the slides into a staining dish of hematoxylin for 30 seconds. Remove and place into an acid bath (200 mL UltraPure water and 1–3 drops of acetic acid). Rinse with UltraPure water.
18. Add several drops of coverslip solution.
19. Place the coverslip on top of the slide.
20. Store slides at room temperature.

---

Rockland products are for research use only and are not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Immunochemicals, Inc., P.O. BOX 5199, Limerick, Pennsylvania 19468, USA.

© 2025 Rockland Immunochemicals, Inc. All rights reserved.