

Datasheet for D405-06-1000**Guinea Pig Plasma Non-Sterile In Sodium Heparin****Overview**

Description:	Guinea Pig Plasma (Non-Sterile) in Sodium Heparin - D405-06-1000
Item No.:	D405-06-1000
Size:	1 L
Origin:	Guinea Pig

Product Details

Synonyms:	guinea pig plasma in Sodium Heparin, non-sterile plasma from guinea pig with anticoagulant Sodium Heparin, non-sterile guinea pig plasma
Species of Origin:	Guinea Pig

Application Details

Application Note:	pH: normal Immunoelectrophoresis: normal Hemoglobin: normal IgG Concentration: normal
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

Tissue Data

Tissue Type:	Plasma
Sex:	Mixed
Strain:	Guinea Pig - Mixed

Formulation

Physical State:	Liquid
Anticoagulant:	Sodium Heparin
Sterility:	Non-sterile
Preservative:	None
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store container at -20° C prior to opening. Avoid cycles of freezing and thawing.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

