

Datasheet for KHG001**Hematoxylin Counterstain Maxtag™ Histo For IHC****Overview**

Description:	Hematoxylin Counterstain MaxTag™ Histo for IHC - KHG001
Item No.:	KHG001
Size:	15 mL
Applications:	IHC

Product Details

Background:	Hematoxylin counterstain reagent is specifically designed for the use in Rockland MaxTag™ Histo for Immunohistochemistry Kit.
Synonyms:	Hematoxylin Counterstain Maxtag™ Histo For IHC, Hematoxylin Counterstain for immunohistochemistry
Detection Kit Type:	Immunohistochemistry Kit

Application Details

Tested Applications:	IHC
Other:	Expiration Date: 03 MAY 2024

Formulation


Physical State:	Liquid
Preservative:	None
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	See kit insert for complete instructions.

Expiration: See kit insert for complete instructions.

Images



Immunohistochemistry

Breast tissue stained with hematoxylin and eosin staining. Used when a stronger or darker nuclear stain is required for cytology or immunohistochemistry (IHC) counterstaining.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

