

**Datasheet for R100-0100****Bovine/Calf Blood Defibrinated****Overview**

<b>Description:</b>	Bovine/Calf Blood Defibrinated - R100-0100
<b>Item No.:</b>	R100-0100
<b>Size:</b>	100 mL
<b>Applications:</b>	LFA, Other
<b>Origin:</b>	Bovine

**Product Details**

<b>Background:</b>	Defibrinated bovine blood is aseptically collected and prepared to order from our colony of animals. Fibrin is mechanically removed during the clotting process without the presence of anticoagulants or other additives. Defibrinated bovine blood is ideal for blood agar plate preparations. Blood Agar is a general purpose enriched medium often used to grow fastidious organisms and to differentiate bacteria based on their hemolytic properties. Defibrinated bovine blood is perishable and is collected and processed upon receipt of your order.
<b>Synonyms:</b>	Sterile Bovine Blood Defibrinated, Sterile Calf Blood Defibrinated, Sterile Bovine/Calf Blood Defibrinated
<b>Species of Origin:</b>	Bovine

**Application Details**

<b>Suggested Applications:</b>	LFA, Other (Based on references)
<b>Application Note:</b>	Blood agar plate preparations
<b>Assay Dilutions:</b>	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

**Tissue Data**

<b>Tissue Type:</b>	Blood
<b>Sex:</b>	Mixed

**Strain:** Bovine/Calf - Mixed

---

## Formulation

**Physical State:** Liquid

---

**Buffer:** None

---

**Anticoagulant:** Defibrinated

---

**Sterility:** Sterile

---

**Preservative:** None

---

**Stabilizer:** None

---

## Shipping & Handling

**Shipping Condition:** Wet Ice

---

**Storage Condition:** Store defibrinated bovine blood at 4° C prior to opening. Use aseptic technique to maintain sterility when opening product. Be advised that blood is a perishable product and exact shelf may depend on application.

---

**Expiration:** This product MAY be stable for up to four (4) weeks if properly stored and handled.

---

## References

- Nileshi Saraf et al. Multiplex Viral Detection Platform Based on a Aptamers-Integrated Microfluidic Channel. *ACS Omega*. (2019)
- Tian Y et al. A disposable blood cyanide sensor. *Anal Chim Acta*. (2013)

## Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

[www.rockland.com](http://www.rockland.com)  
[tech@rockland.com](mailto:tech@rockland.com)  
+1 484.791.3823