

Datasheet for R214-0100**Sterile Human Blood In Sodium Citrate****Overview**

Description:	Sterile Human Blood in Sodium Citrate - Neg for HIV-I antigens, HCV, STS and HBsAg by FDA Licensed tests. Expires xx/xx/xxxx - R214-0100
Item No.:	R214-0100
Size:	100 mL
Origin:	Human

Product Details

Background:	Citrated human blood is aseptically collected. Human red blood cells are perishable and are collected and processed upon receipt of your order. Citrated human blood is ideal for blood agar plate preparations. Blood Agar is a general purpose enriched medium often used to grow fastidious organisms and to differentiate bacteria based on their hemolytic properties. Citrated human RBCs are also useful for the titration of complement, adsorption procedures, testing for agglutinins/HA assays, and for the preparation of stroma as particulate reagents.
Synonyms:	Human Whole Blood, Citrated Human Blood, Human RBCs, Human Red Blood Cells
Species of Origin:	Human

Application Details

Application Note:	Blood agar plates, washed pooled cells preparations, complement titration, adsorption procedures, HA assays and for the preparation of stroma as particulate reagents.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

Tissue Data

Tissue Type:	Blood
Sex:	Mixed
Strain:	Adult

Formulation

Physical State:	Liquid
Buffer:	None
Anticoagulant:	Sodium Citrate
Sterility:	Sterile
Preservative:	None
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Wet Ice
Storage Condition:	Store citrated human blood at 4° C prior to opening. Use aseptic technique to maintain sterility when opening product. Be advised that blood is a perishable product and exact shelf may depend on application.
Expiration:	This product MAY be stable for up to four (4) weeks if properly stored and handled.

Disclaimer

No test method can provide total assurance that the hepatitis B virus, hepatitis C virus, human immunodeficiency virus, or any other infectious agents are absent. Thus, all blood products, including purified proteins derived from human blood sources, should be handled at Biosafety Level 2 as recommended by the CDC\NIH manual entitled Biosafety in Microbiological and Biomedical Laboratories for potentially infectious human serum, blood specimens or proteins derived from same. Source material for the human blood product supplied to your facility has been tested for the detection of HIV antibody, Hepatitis B surface antigen, antibody to Hepatitis C, HIV 1 antigen(s), antibody to HTLV - I/II, and syphilis by FDA guidelines. All units were found to be non-reactive/negative for these tests. All human blood source material is collected in FDA licensed centers and is tested with FDA approved test kits.

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

