

Datasheet for TMS-06

6 Position TrueBlot® Magnetic Separator**Overview**

Description:	6 Position Microcentrifuge Tube TrueBlot® Magnetic Separator - TMS-06
Item No.:	TMS-06
Size:	1 Each

Product Details

Background:	A magnetic separation unit designed to accommodate up to six 1.5 mL microcentrifuge tubes in order to separate magnetic particles contained in those tubes. Optimal working volume is up to 1 mL.
Synonyms:	holder, rack, magnets, magnet

Application Details

Application Note:	Magnetic separator unit can be used to hold, and separate magnetic particles.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.


Formulation

Concentration:	n/a
-----------------------	-----

Shipping & Handling

Shipping Condition:	Ambient
----------------------------	---------

Images

**Separator**

This magnetic separator is designed to accommodate up to six 1.5 mL micro centrifuge tubes.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

