

Datasheet for 001-G70-MF2

Ultrapure Collagen III for plate coating - Bovine Placenta**Overview**

Description:	Ultrapure Collagen III for plate coating - Bovine Placenta - 001-G70-MF2
Item No.:	001-G70-MF2
Size:	8 mL
Applications:	Cellular Assay
Origin:	Bovine

Product Details

Background:	Bovine Collagen Type III is one of a family of proteins found particularly in the flesh and connective tissues of mammals (approximately one-third of the body's total protein). Over two dozen types of collagen have been described, though Type III is the collagen associated with granulation tissues formed during healing.
Synonyms:	Type III collagen, collagen 3, bovine collagen, plate coating solution
Species of Origin:	Bovine

Target Details

Gene Name:	COL3A
Purity/Specificity:	Ultrapure Bovine Collagen Type III Solution is purified from bovine placenta and shown to be 99% pure by gel electrophoresis and immunoassays.
Relevant Links:	Collagen_Top_Culture-Protocol

Application Details

Suggested Applications:	Cellular Assay (Based on references)
Application Note:	Ultrapure Bovine Collagen Type III solution has been tested in cell culture and is suitable for plate coating. Enough to coat 10 (ten) T75 plates.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.


Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1.0mg/mL by dry weight
Buffer:	0.013M HCl
Preservative:	None
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Wet Ice
Storage Condition:	Store vial at 4° C prior to opening. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. Keep under sterile conditions.
Expiration:	Expiration date is six (6) months from date of receipt.

Images



Plate

Ultrapure Collagen III from Bovine Placenta for plate coating. MDA-MB-231 cells grown under cell specific optimal conditions for 24hs on plates coated with 001-G70-MF2 (A) or no coating (B).

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.