

Datasheet for 009-001-B96**IL-9 Recombinant Human Protein****Overview**

Description:	IL-9 Recombinant Human Protein - 009-001-B96
Item No.:	009-001-B96
Size:	10 µg
Applications:	SDS-PAGE
Origin:	Human
Expressed in:	E. coli

Product Details

Background:	Interleukin (IL-9) is known to regulate many cell types involved in T- helper type 2 responses classically associated with asthma, including B- and T- lymphocytes, mast cells, eosinophils and epithelial cells. Growing evidence obtained from human genomic analysis and antigen-challenged transgenic mice suggests that IL-9 is a candidate factor in immunoglobulin E (IgE) production and thus is thought to be associated with bronchial inflammation and bronchial hyperresponsiveness (BHR). Recombinant Human IL-9 produced in E.coli is a single, non-glycosylated single polypeptide chain containing 127 amino acids and having a molecular mass of 14,004 Daltons.
Synonyms:	Interleukin 9 cytokine, IL-9
Species of Origin:	Human
Expressed in:	E. coli
Type:	Recombinant Protein
Low Endotoxin:	Yes

Target Details

Gene Name:	IL9
Purity/Specificity:	Purity was determined to be greater than 98% as determined by analysis by RP-HPLC and by reducing and non-reducing SDS- PAGE, against known standard.
Relevant Links:	• 009-001-B96 SDS

- [NCBI](#)
- [UniProtKB - Q6NZ75](#)

Application Details

Tested Applications:	SDS-PAGE
Application Note:	IL-9 protein has been tested by SDS-PAGE and is suitable as a control for polyclonal or monoclonal anti-IL-9 in immunological assays.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
Other:	Biological Activity: IL-9 is fully biologically active when compared to a standard. The ED50 as determined by the dose-dependent proliferation of human MO7e cells is < 0.5 ng/ml. Endotoxin Level: Measured by LAL is <0.01ng/μg or <0.1EU/μg.


Formulation

Physical State:	Lyophilized
Buffer:	0.01 M Sodium Phosphate, pH 7.5
Preservative:	None
Stabilizer:	None
Reconstitution Volume:	10μl (10-100μl)
Reconstitution Buffer:	Restore with deionized water (or equivalent)

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	Store vial at 4° C prior to restoration. Dilute only prior to immediate use. Maintain sterility. This product DOES NOT contain preservative. DO NOT VORTEX. We recommend adding a carrier protein such as HSA or BSA to 0.1% (i.e. 1.0 mg/mL) . For best results aliquot contents and freeze at -20° C or colder. Avoid cycles of freezing and thawing. Centrifuge vial before each opening to dislodge contents from the cap and to clarify if contents are not clear after standing at room temperature.
Expiration:	Expiration date is six (6) months from date of receipt.

Images



SDS-PAGE

SDS-PAGE shows bands corresponding to IL-9 (1 μ g) in lane 1 (reduced) and lane 3 (non-reduced, arrowhead). Molecular weight estimation was made by comparison to prestained MW markers, lane 2.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.