

Datasheet for 009-001-R40S

Annexin 5 protein-HIS Epitope**Overview**

Description:	Annexin 5 recombinant protein-HIS Epitope - 009-001-R40S
Item No.:	009-001-R40S
Size:	20 µg
Origin:	Human
Expressed in:	E. coli

Product Details

Background:	Annexin 5 is a member of the annexin family of calcium-dependent phospholipid binding proteins. Annexin 5 is a phospholipase A2 and protein kinase C inhibitory protein with calcium channel activity and a potential role in cellular signal transduction, inflammation, growth and differentiation. In vitro studies have shown that Annexin 5 is a cellular protein which plays a role in the inhibition of blood coagulation by competing for phosphatidylserine binding sites with prothrombin (1). Through affinity assays Annexin 5 has been localized on cells that expressed phosphatidylserine on the cell surface, a feature found in apoptosis as well as other forms of cell death (2). Annexin 5 Protein is ideal for investigators involved in Signaling Proteins, Cellular Proteins, Apoptosis/Autophagy, Cancer, and Neurobiology research.
Synonyms:	ANXA5, PP4, ANX5, ENX2, Annexin A5, Anchorin CII, Annexin V, Annexin-5, Calphobindin I, CBP-I, Endonexin II, Lipocortin V, Placental anticoagulant protein 4, PAP-I, Thromboplastin inhibitor, Vascular anticoagulant-alpha, VAC-alpha
Species of Origin:	Human
Expressed in:	E. coli
Type:	Recombinant Protein

Target Details

Gene Name:	ANXA5
Purity/Specificity:	Recombinant full-length human Annexin 5 Protein was expressed by E. coli cells using an N-Terminal his epitope. The purity was determined to be >95% by densitometry.
Relevant Links:	• NCBI - NM_001154

Application Details

Application Note:	Annexin 5 Protein is stored in 50mM sodium phosphate, pH 7.0, 300mM NaCl, 150mM imidazole, 0.1mM PMSF, 2mM DTT, 25% glycerol. Annexin 5 Protein is suitable for use in Western Blot. Expect a band approximately ~ 37kDa on specific lysates or tissues. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
WB:	User Optimized

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	0.2 µg/µL
Buffer:	See application note.

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

