

Datasheet for 009-001-S41S**LOXL2 protein-GST fusion****Overview**

Description:	LOXL2 recombinant protein-GST fusion protein - 009-001-S41S
Item No.:	009-001-S41S
Size:	20 µg
Origin:	Human
Expressed in:	Sf9 cells

Product Details

Background:	LOXL2 is a member of the lysyl oxidase gene family which is essential for the biogenesis of connective tissue. LOXL2 encodes an extracellular copper-dependent amine oxidase that catalyzes the first step in the formation of crosslinks in collagens and elastin. LOXL2 may play an important role in developmental regulation, senescence, tumor suppression, cell growth control, and chemotaxis to each member of the family (1). LOXL2 may mediate fibroblast activation through enzymatic crosslinking of fibrillar collagen and local matrix tension, resulting in activation of TGFb1 signaling (2). LOXL2 Protein is ideal for investigators involved in Signaling Proteins, Lysyl Oxidase Proteins, Cancer, and Metabolic Disorder research.
Synonyms:	LOR2, WS9-14, Lysyl oxidase homolog 2, EC 1.4.3.13, Lysyl oxidase-like protein 2, Lysyl oxidase-related protein 2, Lysyl oxidase-related protein WS9-14
Species of Origin:	Human
Expressed in:	Sf9 cells
Type:	Recombinant Protein

Target Details

Gene Name:	LOXL2
Purity/Specificity:	Recombinant human LOXL2 (519-end) was expressed by baculovirus in Sf9 insect cells using an N-terminal Glutathione-S-Transferase fusion protein. The purity was determined to be >75% by densitometry.
Relevant Links:	• NCBI - BC000594

Application Details

Application Note:	LOXL2 Protein is stored in 50mM Tris-HCl, pH 7.5, 50mM NaCl, 10mM glutathione, 0.1mM EDTA, 0.25mM DTT, 0.1mM PMSF, 25% glycerol. LOXL2 Protein is suitable for use in Western Blot. Expect a band approximately ~65kDa on specific lysates or tissues. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
WB:	User Optimized

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	0.1 µg/µL
Buffer:	See application note.

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

