

Datasheet for 009-001-S46S

MDM4 (1-137) protein-HIS Epitope**Overview**

Description:	MDM4(1-137) recombinant protein-HIS Epitope - 009-001-S46S
Item No.:	009-001-S46S
Size:	20 µg
Origin:	Human
Expressed in:	E. coli

Product Details

Background:	MDM4 is a nuclear protein that contains a p53 binding domain at the N-terminus and a RING finger domain at the C-terminus. MDM4 shows structural similarity to p53-binding protein MDM2 and both proteins bind the p53 tumor suppressor protein and inhibit its activity. However, unlike MDM2, MDM4 does not cause nuclear export or degradation of p53. Instead, MDM4 inhibits p53 activity by binding to the transcriptional activation domain of p53. MDM4 is overexpressed in a variety of human cancers (1). Expression level of MDM4 is significantly higher in chronic lymphocytic leukemia. MDM4 is a specific chemotherapeutic target for treating retinoblastoma (2). MDM4 Protein is ideal for investigators involved in Signaling Proteins, Cell Cycle Proteins, AKT/PKB Pathway, Apoptosis/Autophagy, Cancer, Cell Cycle, and Cellular Stress research.
Synonyms:	DKFZp781B1423, HDMX, MDMX, MGC132766, MRP1, Protein Mdm4, Double minute 4 protein, Mdm2-like p53-binding protein, Protein Mdmx, p53-binding protein Mdm4
Species of Origin:	Human
Expressed in:	E. coli
Type:	Recombinant Protein

Target Details

Gene Name:	MDM4
Purity/Specificity:	Recombinant human MDM4 (1-137) was expressed in E. coli cells using an N-Terminal his epitope. The purity was determined to be >75% by densitometry.
Relevant Links:	<ul style="list-style-type: none">• NCBI - NM_002393

Application Details

Application Note:	MDM4 Protein is stored in 50mM sodium phosphate, pH 7.0, 300mM NaCl, 150mM imidazole, 0.1mM PMSF, 0.25mM DTT, 25% glycerol. MDM4 Protein is suitable for use in Western Blot. Expect a band approximately ~17.3kDa on specific lysates or tissues. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
WB:	User Optimized

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	0.2 µg/µL
Buffer:	See application note.

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

