

Datasheet for 009-001-S48S**MGRN1 protein-GST fusion****Overview**

Description:	MGRN1 recombinant protein-GST fusion protein - 009-001-S48S
Item No.:	009-001-S48S
Size:	20 µg
Origin:	Human
Expressed in:	Sf9 cells

Product Details

Background:	MGRN1 or mahogunin ring finger 1, E3 ubiquitin protein ligase is a C3HC4 RING-containing protein with E3 ubiquitin ligase activity in vitro (1). MGRN1 gene encodes 4 mRNA isoforms that contain a central C3HC4 RING domain as the only recognizable protein motif and regions flanking the RING domain and the N-terminal region of the protein are also conserved between vertebrate and invertebrate genome and MGRN1 is highly expressed in brain (2). MGRN1 has been detected in the adrenal cortex in kidney and its role may relate to the trafficking and/or degradation of the melanocortin 2 receptor. MGRN1 Protein is ideal for investigators involved in Signaling Proteins, Ubiquitin Proteins, Angiogenesis, Apoptosis/Autophagy, Cancer, Cell Cycle, Cellular Stress, Inflammation, Invasion/Metastasis, and Neurobiology research.
Synonyms:	CD41B, GP2B, ITGA2B, KIAA0544, RNF156, E3 ubiquitin-protein ligase MGRN1 EC 6.3.2.-, Mahogunin RING finger protein 1, RING finger protein 156
Species of Origin:	Human
Expressed in:	Sf9 cells
Type:	Recombinant Protein

Target Details

Gene Name:	MGRN1
Purity/Specificity:	Recombinant human MGRN1 (2-end) was expressed by baculovirus in Sf9 insect cells using an N-Terminal Glutathione-S-Transferase fusion protein. The purity was determined to be >90% by densitometry.
Relevant Links:	• NCBI - NM_015246

Application Details

Application Note:	MGRN1 Protein is stored in 50mM Tris-HCl, pH 7.5, 50mM NaCl, 10mM glutathione, 0.1mM EDTA, 0.25mM DTT, 0.1mM PMSF, 25% glycerol. MGRN1 Protein is suitable for use in Western Blot and Kinase Assay. Expect a band approximately ~96kDa on specific lysates or tissues. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
WB:	User Optimized

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	0.1 µg/µL
Buffer:	See application note.

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

