

Datasheet for 009-001-T02S**PP2B protein-GST fusion****Overview**

Description:	PP2B recombinant protein-GST fusion protein - 009-001-T02S
Item No.:	009-001-T02S
Size:	20 µg
Origin:	Human
Expressed in:	Sf9 cells

Product Details

Background:	PP2B is a key regulator of the anaphase-promoting factor which is an excellent substrate for protein phosphatase. PP2B is involved in bovine ocular tissues and is postulated to be involved in the immunologic privilege of the cornea, in retinal signal transduction, and in the toxic effects of immunosuppressants on the eye. PP2B acts as an inhibitory constraint on L-LTP that is relieved by PKA, and that this inhibitory constraint acts as a gate to regulate the synaptic induction of L-LTP (1). PP2B plays a role in the transition from short- to long-term memory and also regulates bone formation through an effect on osteoblast differentiation (2). PP2B Protein is ideal for investigators involved in Signaling Proteins, Cellular Proteins, Cancer, Cardiovascular Disease, Cell Cycle, ERK/MAPK Pathway, Metabolic Disorder, Phosphatases, and Ser/Thr Kinases research.
Synonyms:	PPP2B, PPP3CA, CALN, CCN1, CNA1, CALNA, CALNA1, Serine/threonine-protein phosphatase 2B catalytic subunit alpha isoform, CAM-PRP catalytic subunit, Calmodulin-dependent calcineurin A subunit alpha isoform
Species of Origin:	Human
Expressed in:	Sf9 cells
Type:	Recombinant Protein

Target Details

Gene Name:	PPP3CA
Purity/Specificity:	Recombinant full-length human PP2B was expressed by baculovirus in Sf9 insect cells using an N-Terminal Glutathione-S-Transferase fusion protein. The purity was determined to be >85% by densitometry.

Relevant Links:

- [NCBI - NM_001130691](#)

Application Details

Application Note:	PP2B Protein is stored in 50mM Tris-HCl, pH 7.5, 150mM NaCl, 10mM glutathione, 0.1mM EDTA, 0.25mM DTT, 0.1mM PMSF, 25% glycerol. PP2AB Protein is suitable for use in Western Blot. Expect a band approximately ~83kDa on specific lysates or tissues. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
WB:	User Optimized

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	0.1 µg/µL
Buffer:	See application note.

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

