

Datasheet for 009-001-T06S**PRPF19 protein-GST fusion****Overview**

Description:	PRPF19 recombinant protein-GST fusion protein - 009-001-T06S
Item No.:	009-001-T06S
Size:	20 µg
Origin:	Human
Expressed in:	Sf9 cells

Product Details

Background:	PRPF19 or pre-mRNA-processing factor 19 is the human homolog of yeast Pso4. PRPF19 is a gene essential for cell survival, DNA repair and association with terminal deoxynucleotidyl transferase (1). PRPF19 forms a stable complex with SETMAR (Metnase) which is a SET and transposase fusion protein in humans and plays a positive role in double-strand break (DSB) repair. PRPF19 mediates association of SETMAR with sites of DNA damage (2). SNEV is an orthologue of PRPF19 and is an evolutionarily conserved splicing factor whose oligomerization is necessary for spliceosome assembly. PRPF19 Protein is ideal for investigators involved in Signaling Proteins, Ubiquitin Proteins, Apoptosis/Autophagy, and Cell Cycle Proteins research.
Synonyms:	PSO4, SNEV, PRP19, UBOX4, hPSO4, NMP200, Pre-mRNA-processing factor 19, Nuclear matrix protein 200, PRP19/PSO4 homolog, Senescence evasion factor
Species of Origin:	Human
Expressed in:	Sf9 cells
Type:	Recombinant Protein

Target Details

Gene Name:	PRPF19
Purity/Specificity:	Recombinant human PRPF19 (2-end) was expressed by baculovirus in Sf9 insect cells using an N-Terminal Glutathione-S-Transferase fusion protein. The purity was determined to be >90% by densitometry.
Relevant Links:	• NCBI - NM_014502

Application Details

Application Note:	PRPF19 Protein is stored in 50mM Tris-HCl, pH 7.5, 50mM NaCl, 10mM glutathione, 0.1mM EDTA, 0.25mM DTT, 0.1mM PMSF, 25% glycerol. PRPF19 Protein is suitable for use in Western Blot. Expect a band approximately ~74kDa on specific lysates or tissues. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
WB:	User Optimized

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	0.1 µg/µL
Buffer:	See application note.

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store product at –70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

