

Datasheet for 009-001-T56S

STAT5 protein-HIS Epitope**Overview**

Description:	STAT5 recombinant protein-HIS Epitope - 009-001-T56S
Item No.:	009-001-T56S
Size:	20 µg
Origin:	Human
Expressed in:	Sf9 cells

Product Details

Background:	STAT5 (signal transducers and activators of transcription 5) is a member of the STAT family and this protein is phosphorylated by receptor associated protein kinases in response to cytokines and growth factors (1). The differentiation of T helper (Th) cells is regulated by members of the STAT family of signaling molecules. In BCR-ABL-positive cells, STAT5 is constitutively activated by tyrosine phosphorylation. STAT5 activation results in upregulation of Bcl-XL and increased resistance to induction of apoptosis. STAT5 is involved in expression and growth hormone-mediated sexually dimorphic regulation of cytochrome P450 3A10/lithocholic acid 6 beta-hydroxylase (2). STAT5 Protein is ideal for investigators involved in Signaling Proteins, Transcription Proteins, Apoptosis/Autophagy, Cancer, ERK/MAPK Pathway, Inflammation, and JAK/STAT Pathway research.
Synonyms:	MGF, STAT5A, Signal transducer and activator of transcription 5A
Species of Origin:	Human
Expressed in:	Sf9 cells
Type:	Recombinant Protein

Target Details

Gene Name:	STAT5
Purity/Specificity:	Recombinant full-length human STAT5 was expressed by baculovirus in Sf9 insect cells using an N-Terminal his epitope. The purity was determined to be >90% by densitometry.
Relevant Links:	<ul style="list-style-type: none">• NCBI - NM_003152

Application Details

Application Note:	STAT5 Protein is stored in 50mM sodium phosphate, pH 7.0, 300mM NaCl, 150mM imidazole, 0.1mM PMSF, 0.25mM DTT, 25% glycerol. STAT5 Protein is suitable for use in Western Blot and Kinase Assay. Expect a band approximately ~92kDa on specific lysates or tissues. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
WB:	User Optimized

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	0.2 µg/µL
Buffer:	See application note.

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.



Order online now!

www.rockland.com
tech@rockland.com
+1 484.791.3823