

Datasheet for 009-001-U25S**UBA2 (SAE2) protein-GST fusion****Overview**

Description:	UBA2(SAE2) recombinant protein-GST fusion protein - 009-001-U25S
Item No.:	009-001-U25S
Size:	20 µg
Origin:	Human
Expressed in:	Sf9 cells

Product Details

Background:	UBA2 (SAE2) or ubiquitin-like modifier activating enzyme 2 forms a heterodimer that functions as a SUMO-activating enzyme for the sumoylation of proteins. The posttranslational modification of proteins by the addition of the small protein SUMO, or sumoylation, regulates protein structure and intracellular localization (1). UBA2 functions in SUMO1 activation in a manner analogous to the single E1 ubiquitin-activating enzymes (2). The inactivation of UBA2 may be a therapeutic strategy in MYC-driven cancers. UBA2 Protein is ideal for investigators involved in Signaling Proteins, Ubiquitin Proteins, Cancer, Cell Cycle, and Neurobiology research.
Synonyms:	UBA2, UBLE1B, SAE2, ARX, HRIHFB2115, SUMO-activating enzyme subunit 2
Species of Origin:	Human
Expressed in:	Sf9 cells
Type:	Recombinant Protein

Target Details

Gene Name:	UBA2
Purity/Specificity:	Recombinant full-length human UBA2 (SAE2) was expressed by baculovirus in Sf9 insect cells using an N-Terminal Glutathione-S-Transferase fusion protein. The purity was determined to be >85% by densitometry.
Relevant Links:	NCBI - NM_005499

Application Details

Application Note:	UBA2 Protein is stored in 50mM Tris-HCl, pH 7.5, 50mM NaCl, 10mM glutathione, 0.1mM EDTA, 0.25mM DTT, 0.1mM PMSF, 25% glycerol. UBA2 Protein is suitable for use in Western Blot. Expect a band approximately ~115kDa on specific lysates or tissues. Specific conditions for reactivity should be optimized by the end user.
--------------------------	---

Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
-------------------------	---

WB:	User Optimized
------------	----------------

Formulation

Physical State:	Liquid (sterile filtered)
------------------------	---------------------------

Concentration:	0.1 µg/µL
-----------------------	-----------

Buffer:	See application note.
----------------	-----------------------

Shipping & Handling

Shipping Condition:	Dry Ice
----------------------------	---------

Storage Condition:	Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.
---------------------------	---

Expiration:	Expiration date is one (1) year from date of receipt.
--------------------	---

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

