

Datasheet for 009-001-U29S**UBE2C protein-HIS Epitope****Overview**

Description:	UBE2C recombinant protein-HIS Epitope - 009-001-U29S
Item No.:	009-001-U29S
Size:	20 µg
Reactivity:	H. sapiens (Human)
Origin:	Human
Expressed in:	E. coli

Product Details

Background:	UBE2C or ubiquitin-conjugating enzyme E2C is a member of the E2 ubiquitin-conjugating enzyme family which is required for the destruction of mitotic cyclins and for cell cycle progression (1). UBE2C acted as a dominant-negative inhibitor which blocks the destruction of mitotic cyclins A and B and the onset of anaphase in frog embryos and mammalian cells (2). High UBE2C mRNA expression is associated with poor disease-free survival and overall survival. High tumor grade, as well as high Ki67 protein expression, was more frequent in the high-expression group of UBE2C. UBE2C Protein is ideal for investigators involved in Signaling Proteins, Ubiquitin Proteins, Cancer, Cell Cycle, and Neurobiology research.
Synonyms:	dJ447F3.2, UBCH10, Ubiquitin-conjugating enzyme E2 C, Ubiquitin carrier protein C, Ubiquitin-protein ligase C
Species of Origin:	Human
Expressed in:	E. coli
Type:	Recombinant Protein

Target Details

Gene Name:	UBE2C
Reactivity:	H. sapiens (Human)
Purity/Specificity:	Recombinant full length human UBE2C was expressed in E. coli cells using an N-Terminal his epitope. The purity was determined to be >95% by densitometry.
Relevant Links:	• NCBI - NM_007019

Application Details

Application Note:	UBE2C Protein is stored in 50mM sodium phosphate, pH 7.0, 300mM NaCl, 150mM imidazole, 0.1mM PMSF, 0.25mM DTT, 25% glycerol. UBE2C Protein is suitable for use in Western Blot. Expect a band approximately ~21kDa on specific lysates or tissues. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
WB:	User Optimized

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	0.1 µg/µL
Buffer:	See application note.

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

