

Datasheet for 009-001-U53S**UHKM1 (KIS) protein-GST fusion****Overview**

Description:	UHKM1 (KIS) recombinant protein-GST fusion protein - 009-001-U53S
Item No.:	009-001-U53S
Size:	20 µg
Origin:	Human
Expressed in:	Sf9 cells

Product Details

Background:	UHKM1 (KIS) or U2AF homology motif (UHM) kinase 1 is a serine/threonine protein kinase that promotes cell cycle progression through G1 phase. UHKM1 phosphorylates the cyclin-dependent kinase inhibitor 1B (p27Kip1) and this results in its export to the nucleus where it undergoes degradation. UHKM1 is a growth factor-dependent nuclear kinase which can regulate cell cycle progression (1). Elevated levels of UHKM1 protein in leukemia cells has been shown to promote cell cycle progression. UHKM1 has also been shown to function in the adult nervous system where it is highly expressed and the gene has been associated with schizophrenia (2). UHKM1 Protein is ideal for investigators involved in Signaling Proteins, Cellular Proteins, Cancer, Cell Cycle, and Ser/Thr Kinases research.
Synonyms:	UHKM1, KIS, KIST, Serine/threonine-protein kinase Kist, Kinase interacting with stathmin, PAM COOH-terminal interactor protein 2, P-CIP2, U2AF homology motif kinase 1
Species of Origin:	Human
Expressed in:	Sf9 cells
Type:	Recombinant Protein

Target Details

Gene Name:	UHKM1
Purity/Specificity:	Recombinant full-length human UHKM1 was expressed by baculovirus in Sf9 insect cells using an N-terminal Glutathione-S-Transferase fusion protein. The purity was determined to be >70% by densitometry.
Relevant Links:	• NCBI - NM_144624

Application Details

Application Note:	UHK1 Protein is stored in 50mM Tris-HCl, pH 7.5, 50mM NaCl, 10mM glutathione, 0.1mM EDTA, 0.25mM DTT, 0.1mM PMSF, 25% glycerol. UHK1 Protein is suitable for use in Western Blot. Expect a band approximately ~65kDa on specific lysates or tissues. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
WB:	User Optimized

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	0.05 µg/µl
Buffer:	See application note.

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

