

Datasheet for 009-001-U64-0005**rHuman Artemin Protein****Overview**

Description:	Human Artemin Recombinant Protein - 009-001-U64-0005
Item No.:	009-001-U64-0005
Size:	5 µg
Applications:	SDS-PAGE
Origin:	Human
Expressed in:	E. coli

Product Details

Background:	Artemin is a novel member of the glial cell line-derived neurotrophic factor (GDNF) ligand family. Current evidence suggests that Artemin signals through the receptor complex, GFR α 3-RET, to influence neuron survival in vitro and in vivo. Recombinant human Artemin is a non-glycosylated, disulfide-linked homodimer, containing two 113 amino acid chains, with a total molecular weight of 24.2 kDa.
Synonyms:	ARTN, enovin, neublastin
Species of Origin:	Human
Expressed in:	E. coli
Type:	Recombinant Protein
Low Endotoxin:	Yes

Target Details

Gene Name:	ARTN
Purity/Specificity:	Artemin purity was determined to be greater than 98% as determined by analysis by HPLC and by reducing and non-reducing SDS-PAGE.
Relevant Links:	• UniProtKB - Q5T4W7

Application Details

Tested Applications:	SDS-PAGE
Application Note:	Recombinant human Artemin has been tested by SDS-PAGE and is suitable as a control for polyclonal or monoclonal anti-Artemin in immunological assays.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
Other:	Endotoxin Level: Measured by kinetic LAL analysis and is typically ≤ 1 EU/ μ g protein. Biologic Activity: The activity is determined by the dose-dependant proliferation of the SH-SY5Y cell line and is typically 4-8 ng/mL. The activity can also be determined by its ability to promote survival and neurite outgrowth.


Formulation

Physical State:	Lyophilized
Buffer:	0.1% Trifluoroacetic acid
Preservative:	None
Stabilizer:	None
Reconstitution Volume:	5 μ l (5-50 μ l)
Reconstitution Buffer:	Restore with deionized water (or equivalent)

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	Store vial at 4° C prior to restoration. Dilute only prior to immediate use. Maintain sterility. This product DOES NOT contain preservative. DO NOT VORTEX. We recommend adding a carrier protein such as HSA or BSA to 0.1% (i.e. 1.0 mg/mL). For best results aliquot contents and freeze at -20° C or colder. Avoid cycles of freezing and thawing. Centrifuge vial before each opening to dislodge contents from the cap and to clarify if contents are not clear after standing at room temperature.
Expiration:	Expiration date is six (6) months from date of receipt.

Images



SDS-PAGE

SDS-PAGE of Human Artemin Recombinant Protein. Lane 1: 1 μ g Human Artemin in non-reducing conditions (-). Lane 2: Molecular weight marker. Lane 3: 1 μ g Human Artemin in reducing conditions (+). Lane 4: Molecular weight marker. Human Artemin is a predicted homodimer with a MW of 24 kDa.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.