

**Datasheet for 009-001-U74-0020****rHuman Endocrine Gland VEGF Protein****Overview**

|                      |                                                                                                 |
|----------------------|-------------------------------------------------------------------------------------------------|
| <b>Description:</b>  | Human Endocrine Gland Vascular Endothelial Growth Factor Recombinant Protein - 009-001-U74-0020 |
| <b>Item No.:</b>     | 009-001-U74-0020                                                                                |
| <b>Size:</b>         | 20 µg                                                                                           |
| <b>Applications:</b> | SDS-PAGE                                                                                        |
| <b>Origin:</b>       | Human                                                                                           |
| <b>Expressed in:</b> | E. coli                                                                                         |

**Product Details**

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Background:</b>        | Endocrine Gland-derived Vascular Endothelial Growth Factor (EG-VEGF) is an angiogenic growth factor specifically expressed in the ovaries, testis, adrenal and placental tissues. The identification of tissue-selective angiogenic factors raises the possibility that other secreted molecules in this class exist. EG-VEGF expression correlates with vascularity in polycystic ovary syndrome, a leading cause of infertility. Recombinant human EG-VEGF is a non-glycosylated protein, containing 86 amino acids, with a molecular weight of 9.6 kDa. |
| <b>Synonyms:</b>          | Prokineticin 1,                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Species of Origin:</b> | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Expressed in:</b>      | E. coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Type:</b>              | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Low Endotoxin:</b>     | Yes                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |

**Target Details**

|                            |                                                                                                                                                                                                        |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name:</b>          | PROK1                                                                                                                                                                                                  |
| <b>Purity/Specificity:</b> | Endocrine Gland Vascular Endothelial Growth Factor purity was determined to be greater than 97% as determined by analysis by HPLC, UV-Spectroscopy at 280nm and by reducing and non-reducing SDS-pAGE. |
| <b>Relevant Links:</b>     | <ul style="list-style-type: none"><li>• <a href="#">UniProtKB - P58294</a></li></ul>                                                                                                                   |

## Application Details

|                             |                                                                                                                                                                                                                                               |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Tested Applications:</b> | SDS-PAGE                                                                                                                                                                                                                                      |
| <b>Application Note:</b>    | Endocrine Gland Vascular Endothelial Growth Factor Recombinant Protein has been tested by SDS-PAGE and is suitable as a control for polyclonal or monoclonal anti-Endocrine Gland Vascular Endothelial Growth Factor in immunological assays. |
| <b>Assay Dilutions:</b>     | All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.                                                                                                                                               |
| <b>Other:</b>               | Endotoxin Level: Measured by kinetic LAL analysis and is typically $\leq 1$ EU/ $\mu$ g protein. Biologic Activity: The activity is determined by the dose-dependent proliferation of MIA PaCa-2 cells and is typically 1-4 $\mu$ g/mL.       |


## Formulation

|                               |                                              |
|-------------------------------|----------------------------------------------|
| <b>Physical State:</b>        | Lyophilized                                  |
| <b>Concentration:</b>         | 0.1mg/ml                                     |
| <b>Buffer:</b>                | 0.1% Trifluoroacetic acid                    |
| <b>Preservative:</b>          | None                                         |
| <b>Stabilizer:</b>            | None                                         |
| <b>Reconstitution Volume:</b> | 20 $\mu$ l (20-200 $\mu$ l)                  |
| <b>Reconstitution Buffer:</b> | Restore with deionized water (or equivalent) |

## Shipping & Handling

|                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Shipping Condition:</b> | Ambient                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Storage Condition:</b>  | Store vial at 4° C prior to restoration. Dilute only prior to immediate use. Maintain sterility. This product DOES NOT contain preservative. DO NOT VORTEX. We recommend adding a carrier protein such as HSA or BSA to 0.1% (i.e. 1.0 mg/mL). For best results aliquot contents and freeze at -20° C or colder. Avoid cycles of freezing and thawing. Centrifuge vial before each opening to dislodge contents from the cap and to clarify if contents are not clear after standing at room temperature. |
| <b>Expiration:</b>         | Expiration date is six (6) months from date of receipt.                                                                                                                                                                                                                                                                                                                                                                                                                                                   |

## Images

**SDS-PAGE**

SDS-PAGE of Human Endocrine Gland Vascular Endothelial Growth Factor Recombinant Protein. Lane 1: 1 µg Human EG-VEGF in reducing conditions. Molecular weight marker. Human EG-VEGF is predicted to have a MW of 9.6 kDa.

**Disclaimer**

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.