

Datasheet for 010-001-S39S**LMTK3 protein GST fusion****Overview**

Description:	LMTK3 recombinant protein GST fusion protein - 010-001-S39S
Item No.:	010-001-S39S
Size:	20 µg
Origin:	Mouse
Expressed in:	Sf9 cells

Product Details

Background:	LMTK3, or Lemur tyrosine kinase 3, belongs to the LMR subfamily of tyr protein kinases. LMTK3 decreases the activity of protein kinase C (PKC) and also leads to phosphorylation of AKT (Ser473), thereby increasing binding of forkhead box O3 (FOXO3) to the ESR1 promoter. LMTK3 phosphorylates ERα, protecting it from proteasomal degradation in vitro. Knockdown of LMTK3 leads to a reduction of tumor volume in an orthotopic mouse model and abrogates proliferation of ERα+ but not ERα- cells, indicative of its role in ERα activity. These findings reveal LMTK3 as a new therapeutic target for breast cancer (1). LMTK3 Protein is ideal for investigators involved in Signaling Proteins, Cellular Proteins, Cancer, and Ser/Thr Kinases research.
Synonyms:	AATYK3, MGC90743, Serine/threonine-protein kinase LMTK3, EC 2.7.11.1, Apoptosis-associated tyrosine kinase 3, Lemur tyrosine kinase 3
Species of Origin:	Mouse
Expressed in:	Sf9 cells
Type:	Recombinant Protein

Target Details

Gene Name:	Lmtk3
Purity/Specificity:	Recombinant mouse LMTK3 (63-484) was expressed by baculovirus in Sf9 insect cells using an N-terminal Glutathione-S-Transferase fusion protein. The purity was determined to be >75% by densitometry.
Relevant Links:	• NCBI - NM_001005511

Application Details

Application Note:	LMTK3 Protein is stored in 50mM Tris-HCl, pH 7.5, 50mM NaCl, 10mM glutathione, 0.1mM EDTA, 0.25mM DTT, 0.1mM PMSF, 25% glycerol. LMTK3 Protein is suitable for use in Western Blot. Expect a band approximately ~75kDa on specific lysates or tissues. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
WB:	User Optimized

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	0.1 µg/µL
Buffer:	See application note.

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

