

Datasheet for 025-0244**Llama IgG2 isotype control Fluorescein****Overview**

Description:	Llama IgG2 Isotype Control Fluorescein Conjugated - 025-0244
Item No.:	025-0244
Size:	100 µg
Origin:	Llama

Product Details

Background:	Comparative studies of old world camelids (<i>Camelus bactrianus</i> and <i>Camelus dromedarius</i>) and new world camelids (<i>Lama pacos</i>, <i>Lama glama</i> and <i>Lama vicugna</i>) have shown that heavy-chain-only immunoglobulins represent between 35% - 70% of total IgG in the sera of all species. Such antibodies are homodimers of heavy chains that lack the CH1 domain of conventional antibodies and therefore do not interact with light chains, exhibiting a lower molecular weight ~100 kDa. In llama and other species of camelids, these heavy-chain-only immunoglobulins belong to the IgG2 and IgG3 subclasses of gamma chain antibodies.
Synonyms:	Llama IgG2 isotype, Llama IgG2 subclass isotype, Llama IgG2a, IgG2b, IgG2c, Llama IgG2 Fluorescein Conjugated, Llama IgG2 FITC, Llama isotype control
Species of Origin:	Llama
Conjugate:	Fluorescein (FITC)
Format:	IgG
Type:	Native Protein

Target Details

Purity/Specificity:	Llama IgG2 isotype control (which consists of IgG2a, IgG2b and IgG2c isotypes) has been prepared from llama serum by multiple chromatography steps using a combination of protein A and protein G chromatography. Coomassie stained SDS-PAGE of non-reduced unconjugated llama IgG2 shows a band of ~100 kDa whereas the reduced form exhibits ~46 kDa. No bands corresponding to unconjugated llama IgG1 or IgG3 are observed. Fluorescein Conjugated Llama IgG2 was proven by Dot Blot.
----------------------------	--

Application Details

Application Note:	Llama IgG2 FITC conjugated isotype control can be utilized in Flow Cytometry, immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting experiments where determination of sample isotype is important.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
FC:	User Optimized
FLISA:	User Optimized
IF:	User Optimized
WB:	User Optimized

Formulation

Physical State:	Lyophilized
Concentration:	1.0 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Sodium Azide
Stabilizer:	10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free
Reconstitution Volume:	100 µL
Reconstitution Buffer:	Restore with deionized water (or equivalent)

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.