

**Datasheet for 100-401-F38****GRP78 Antibody****Overview**

|                      |                                                                                       |
|----------------------|---------------------------------------------------------------------------------------|
| <b>Description:</b>  | Anti-GRP78 (RABBIT) Antibody - 100-401-F38                                            |
| <b>Item No.:</b>     | 100-401-F38                                                                           |
| <b>Size:</b>         | 200 µL                                                                                |
| <b>Applications:</b> | IF, WB                                                                                |
| <b>Reactivity:</b>   | Human, Mouse, Rat, Bovine, Cryphonectria parasitica, Hamster, Monkey, Rabbit, Xenopus |
| <b>Host Species:</b> | Rabbit                                                                                |

**Product Details**

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Background:</b>   | GRP78 is a ubiquitously expressed, 78-kDa glucose regulated protein, and is commonly referred to as an immunoglobulin chain binding protein (BiP). The BiP proteins are categorized as stress response proteins because they play an important role in the proper folding and assembly of nascent protein and in the scavenging of misfolded proteins in the endoplasmic reticulum lumen. Translation of BiP is directed by an internal ribosomal entry site (IRES) in the 5' nontranslated region of the BiP mRNA. BiP IRES activity increases when cells are heat stressed. GRP78 is also critical for maintenance of cell homeostasis and the prevention of apoptosis. Luo et al. have provided findings that suggest GRP78 is essential for embryonic cell growth and pluripotent cell survival. |
| <b>Synonyms:</b>     | BiP, Grp78, HspA5, MIF2, immunoglobulin heavy chain binding protein, 78 kDa glucose-regulated protein, Heat shock 70 kDa protein 5, Endoplasmic reticulum luminal Ca(2+)-binding protein grp78                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Host Species:</b> | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Clonality:</b>    | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Format:</b>       | Antiserum                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |

**Target Details**

|                        |                                                                                       |
|------------------------|---------------------------------------------------------------------------------------|
| <b>Gene Name:</b>      | Hspa5                                                                                 |
| <b>Reactivity:</b>     | Human, Mouse, Rat, Bovine, Cryphonectria parasitica, Hamster, Monkey, Rabbit, Xenopus |
| <b>Immunogen Type:</b> | Conjugated Peptide                                                                    |

|                            |                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen:</b>          | Grp78 Antibody was produced from whole rabbit serum prepared by repeated immunizations with rat GRP78 (Bip) synthetic peptide conjugated to KLH.                                                                                                                                                                                                                                                                                  |
| <b>Purity/Specificity:</b> | Anti-GRP78 Antibody was prepared from monospecific antiserum by delipidation and defibrination. A BLAST analysis was used to suggest cross-reactivity with GRP78 from Human, Mouse, Rat, Cow, Fungus ( <i>Cryptosporidium parvum</i> ), Hamster, Monkey, Rabbit, <i>Xenopus</i> based on 100% homology with the immunizing sequence. Cross-reactivity with GRP78 from other sources has not been determined. Heat Shock research. |
| <b>Relevant Links:</b>     | <ul style="list-style-type: none"><li>• <a href="#">NCBI - NP_037215.1</a></li><li>• <a href="#">GeneID - 25617</a></li><li>• <a href="#">UniProtKB - P06761</a></li></ul>                                                                                                                                                                                                                                                        |

## Application Details

|                             |                                                                                                                                                                                                                                                   |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Tested Applications:</b> | IF, WB                                                                                                                                                                                                                                            |
| <b>Application Note:</b>    | Anti-GRP78 Antibody is tested in WB and IF and suitable for use in IHC and IP. Expect a band approximately ~78kDa corresponding to the molecular weight of Grp78 protein. Specific conditions for reactivity should be optimized by the end user. |
| <b>Assay Dilutions:</b>     | All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.                                                                                                                                                   |
| <b>IF:</b>                  | 1:100                                                                                                                                                                                                                                             |
| <b>WB:</b>                  | 1:1000                                                                                                                                                                                                                                            |

## Formulation

|                        |                           |
|------------------------|---------------------------|
| <b>Physical State:</b> | Liquid (sterile filtered) |
| <b>Concentration:</b>  | 85 mg/mL by Refractometry |
| <b>Buffer:</b>         | 1X PBS, pH 7.4            |
| <b>Preservative:</b>   | 0.09% (w/v) Sodium Azide  |
| <b>Stabilizer:</b>     | 50% (v/v) Glycerol        |


## Shipping & Handling

|                            |         |
|----------------------------|---------|
| <b>Shipping Condition:</b> | Dry Ice |
|----------------------------|---------|

**Storage Condition:** Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.


**Expiration:** Expiration date is one (1) year from date of receipt.

## Images



### Western Blot

Western Blot of rabbit anti-Grp78 antibody. Lane 1: GRP78 recombinant protein. Lane 2: Human tissues. Lane 3: Dog tissues. Lane 4: Mouse tissues. Load: 10 µg per lane. Primary antibody: Grp78 antibody at 1:1000 for overnight at 4°C. Secondary antibody: IRDye800™ rabbit secondary antibody at 1:10,000 for 45 min at RT. Block: 5% BLOTTO overnight at 4°C. Predicted/Observed size: 72.4 kDa, ~75 kDa for Grp78. Other band(s): none.



### Western Blot

Western Blot of rabbit anti-Grp78 antibody. Lane 1-6: Rat tissues. Load: 10 µg per lane. Primary antibody: Grp78 antibody at 1:1000 for overnight at 4°C. Secondary antibody: IRDye800™ rabbit secondary antibody at 1:10,000 for 45 min at RT. Block: 5% BLOTTO overnight at 4°C. Predicted/Observed size: 72.4 kDa, ~80 kDa for Grp78. Other band(s): none.

## Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.