

Datasheet for 100-401-V89**Wdr5 Antibody****Overview**

Description:	Anti-Wdr5 (RABBIT) Antibody - 100-401-V89
Item No.:	100-401-V89
Size:	100 µL
Applications:	ELISA, WB
Reactivity:	Mouse
Host Species:	Rabbit

Product Details

Background:	Wdr5 belongs to the family of WD repeat proteins which are involved in different cellular processes such as cell cycle progression, signal transduction, apoptosis, and gene regulation. It is a component of the Set1A and Set1B histone H3 methylation complexes. These complexes methylate lysine 4 of H3, thereby activating gene transcription. Wdr5 interacts with H3 dimethyl K4, but not with tri or mono methylated H3K4. Anti-Wdr5 Antibody is ideal for research in Cell Cycle Control, Signal Transduction, Apoptosis, Epigenetics, Gene Expression and Cancer.
Synonyms:	WD repeat-containing protein 5, BMP2-induced 3-kb gene protein, WD repeat-containing protein BIG-3, Big, Big3
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	Antiserum

Target Details

Gene Name:	Wdr5
Reactivity:	Mouse
Immunogen Type:	Conjugated Peptide
Immunogen:	Anti-Wdr5 Antibody was produced in rabbits by repeated immunizations with two synthetic peptides containing an amino acid sequence from the internal part of mouse Wdr5.

Purity/Specificity: Anti-Wdr5 Antibody is monospecific antiserum processed by delipidation and defibrination followed by sterile filtration. Cross reactivity with Wdr5 from other species not tested.

Relevant Links:

- [UniProtKB - P61965](#)
- [GeneID - 140858](#)
- [NCBI - NP_543124.1](#)

Application Details

Tested Applications: ELISA, WB

Application Note: Anti-Wdr5 Antibody has been tested by ELISA and Western Blots. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 36 kDa in the appropriate cell lysate or extract.

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

ELISA: 1:100 - 1:500

WB: 1:1,000

Formulation

Physical State: Liquid (sterile filtered)

Concentration: neat

Buffer: None

Preservative: 0.05% (w/v) Sodium Azide

Stabilizer: None



Shipping & Handling

Shipping Condition: Dry Ice

Storage Condition: Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Expiration: Expiration date is one (1) year from date of receipt.

Images



ELISA

ELISA results of Rabbit anti-Wdr5 antibody. Antigen: BSA conjugated Wdr5. Coating amount: 0.1 µg per well. Dilution series: serial dilution. Estimated Antibody Titer to be 1:12,300. Substrate: TMB (p/n TMBE-1000).

Western Blot

Western Blot results of Rabbit anti-Wdr5 antibody. Lane 1: NIH3T3 nuclear extracts lysate. Lane 2: embryonic stem cells E14Tg2a lysate. Load: 20 µg per lane. Primary antibody: Wdr5 antibody at 1:1000 for overnight at 4°C. Secondary antibody: goat anti-rabbit HRP secondary antibody at 1:10,000 for 45 min at RT. Block: BSA/PBS-Tween overnight at 4°C. Predicted: 36kDa.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.