

Datasheet for 100-4159

Hexokinase Antibody

Overview

Description:	Anti-Hexokinase (Yeast) (RABBIT) Antibody - 100-4159
Item No.:	100-4159
Size:	2 mL
Applications:	WB
Reactivity:	Hexokinase (Yeast)
Host Species:	Rabbit

Product Details

Background:	Hexokinase is the main glucose phosphorylating enzyme. It may play a regulatory role in both induction and repression of gene expression by glucose. In yeast there are three glucose-phosphorylating isoenzymes, designated hexokinase I, II and glucokinase. This protein is involved in the pathway hexose metabolism, which is part of Carbohydrate metabolism.
Synonyms:	rabbit anti-Hexokinase Antibody, DKFZp686M1669 antibody, Hexokinase 2 antibody, Hexokinase 2 muscle antibody, Hexokinase type II antibody
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	Antiserum

Target Details

Gene Name:	HXK2
Reactivity:	Hexokinase (Yeast)
Immunogen Type:	Native Protein
Immunogen:	Hexokinase [Yeast]
Purity/Specificity:	This product was prepared from monospecific antiserum by a delipidation and defibrination. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-rabbit serum, purified and partially purified Hexokinase [Yeast]. Cross reactivity against Hexokinase from other tissues and species may occur but have not been specifically determined.

Relevant Links:	• UniProtKB - P04806 • NCBI - CAA96973.1 • UniProtKB - P04807 • GeneID - 852639
------------------------	---

Application Details

Suggested Applications:	WB (Based on references)
Application Note:	Anti-Hexokinase is suitable for use in ELISA, western blot, and immunohistochemistry. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:5,000 - 1:25,000
IHC:	User Optimized
WB:	1:500 - 1:3,000


Formulation

Physical State:	Lyophilized
Concentration:	80 mg/mL by Refractometry
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Sodium Azide
Stabilizer:	None
Reconstitution Volume:	2.0 mL
Reconstitution Buffer:	Restore with deionized water (or equivalent)

Shipping & Handling


Shipping Condition:	Ambient
Storage Condition:	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images




Western Blot

Repression of HML and HMR preceded heterochromatin maturation. (B) Protein immunoblotting in a strain expressing sir3-8-3xV5 (JRY13467) constitutively at 25°C (first lane), constitutively at 37°C (second lane), and at 30, 60, 90, 120, and 150 min after a shift to 25°C. Top row is 3xV5-tagged sir3-8 protein, the middle row is the same as the top row but at a higher exposure, and the bottom row is the loading control Hxk2. Figure provided by CiteAb. Source: Elife, PMID: 35073254.



Western Blot

Nucleosome binding was required for spread, but not recruitment, of Sir3 to regions of heterochromatin. (B) Protein immunoblotting in strains expressing Sir3 (no tag, JRY11699), Sir3-3xV5 (JRY12601), and sir3-bahΔ-3xV5 (JRY13621). Top row are 3xV5-tagged Sir3 proteins, and bottom row is the loading control Hxk2. Figure provided by CiteAb. Source: Elife, PMID: 35073254.



Western Blot

Effects of cohesin depletion and tT(AGU)C deletion on silencing establishment. (A) SIR3-EBD strains (JRY12269, JRY12270) and SIR3-EBD strains with seamless deletion of tT (AGU)C (JRY12267; JRY12268) were arrested in G1 with α factor, then split, with half the culture receiving estradiol and half receiving ethanol. Samples were collected after 3 hr for RT-qPCR, with each sample normalized to its own pre-estradiol value. (B) Cells with MCD1-AID (JRY12560, JRY12561) were arrested in G1 with α factor, then split, with half receiving auxin and the other half receiving DMSO (solvent control). After 30 min, each culture was further split, with half receiving estradiol and the other half ethanol. All cultures were released to G2/M by addition of protease and nocodazole. Cells were collected after 3 hr for RT-qPCR, with each sample normalized to its own pre-estradiol value. (C) Immunoblot analysis showing Mcd1-AID depletion for experiment described in (B). Figure provided by CiteAb. Source: Elife, PMID: 32687055.

References

- Ben-Menachem RH et al. Mitochondrial-derived vesicles retain membrane potential and contain a functional ATP synthase. *EMBO Rep.* (2023)
- Saxton DS et al. Distinct silencer states generate epigenetic states of heterochromatin. *Mol Cell.* (2022)
- Brothers, M et al. Distinguishing between recruitment and spread of silent chromatin structures in *Saccharomyces cerevisiae*. *ELife* (2022)
- Otto GM et al. Programmed cortical ER collapse drives selective ER degradation and inheritance in yeast meiosis. *J Cell Biol.* (2021)
- Farris D et al. A novel allele of SIR2 reveals a heritable intermediate state of gene silencing. *Genetics.* (2021)
- Goodnight D et al. S-phase-independent silencing establishment in *Saccharomyces cerevisiae*. *Elife.* (2020)
- Eisenberg AR et al. Translation initiation site profiling reveals widespread synthesis of non-AUG-initiated protein isoforms in yeast. *Cell Syst.* (2020)
- Muller R et al. CiBER-seq dissects genetic networks by quantitative CRISPRi profiling of expression phenotypes. *Science.* (2020)
- Lim G et al. Phosphoregulation of Rad51/Rad52 by CDK1 functions as a molecular switch for cell cycle–specific activation of homologous recombination. *Sci Adv.* (2020)
- Sawyer EM et al. Developmental regulation of an organelle tether coordinates mitochondrial remodeling in meiosis. *J Cell Biol.* (2019)
- King GA et al. Meiotic cellular rejuvenation is coupled to nuclear remodeling in budding yeast. *Elife.* (2019)
- Brothers M et al. Mutations in the PCNA DNA Polymerase Clamp of *Saccharomyces cerevisiae* Reveal Complexities of the Cell Cycle and Ploidy on Heterochromatin Assembly. *Genetics.* (2019)
- Brace JL et al. A cell separation checkpoint that enforces the proper order of late cytokinetic events. *J Cell Biol.* (2019)
- Lam, MHY et al. A Comprehensive Membrane Interactome Mapping of Sho1p Reveals Fps1p as a Novel Key Player in the Regulation of the HOG Pathway in *S. cerevisiae*. *Journal of Molecular Biology* (2015)
- Thanabalu T et al. Verprolin function in endocytosis and actin organization: Roles of the Las17p (yeast WASP)-binding domain and a novel C-terminal actin-binding domain. *FEBS J.* (2007)

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.