

Datasheet for 200-101-230S

Maltose Phosphorylase Antibody**Overview**

Description:	Anti-Maltose Phosphorylase (E. coli) (GOAT) Antibody - 200-101-230S
Item No.:	200-101-230S
Size:	25 µL
Applications:	WB
Reactivity:	Bacteria
Host Species:	Goat

Product Details

Background:	Maltose Phosphorylase catalyzes the phosphorolysis of maltose, leading to the formation of glucose and glucose 1-P.
Synonyms:	goat anti-Maltose Phosphorylase Antibody, mapA antibody, mpA antibody
Host Species:	Goat
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	LAF
Reactivity:	Bacteria
Immunogen Type:	Native Protein
Immunogen:	Maltose Phosphorylase [E.coli]
Purity/Specificity:	Anti-MALTOSE PHOSPHORYLASE is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum as well as purified and partially purified Maltose Phosphorylase [E.coli]. Cross reactivity against Maltose Phosphorylase from other sources is unknown.

Relevant Links:

- [UniProtKB - B2GEX2](#)
- [NCBI - YP_001842841.1](#)
- [GeneID - 6233598](#)

Application Details

Tested Applications:	WB
Application Note:	Anti-Maltose Phosphorylase has been tested by western blot. A working dilution of 1:2,000 to 1:10,000 of the reconstitution concentration is suggested for this product.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:5,000 - 1:20,000
IP:	1:100
WB:	1:500 - 1:5,000


Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Sodium Azide
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store vial at -20° C or below prior to opening. This vial contains a relatively low volume of reagent (25 µL). To minimize loss of volume dilute 1:10 by adding 225 µL of the buffer stated above directly to the vial. Recap, mix thoroughly and briefly centrifuge to collect the volume at the bottom of the vial. Use this intermediate dilution when calculating final dilutions as recommended below. Store the vial at -20°C or below after dilution. Avoid cycles of freezing and thawing.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



Western Blot

Western Blot Results of Goat Anti-Maltose Phosphorylase Antibody.

Lane 1: Maltose Phosphorylase Reduced [0.1µg].

Lane 2: Opal Prestained Molecular Weight [p/n MB-210-0500].

Lane 3: Maltose Phosphorylase Non-Reduced [0.1µg].


Primary Antibody: Anti-Maltose Phosphorylase [Goat]
Antibody at 1:1000 overnight at 2-8°C.

Secondary Antibody: Anti-Goat IgG [Rabbit] Peroxidase
Conjugated [p/n 605-4302] at 1:40,000 for 30mins at RT.

Block: Blocking Buffer for Fluorescent Western Blotting [p/n MB-070] for 1hr at RT.

Expect: ~90kDa.

Exposure: 3.4sec.



Western Blot

Rockland Goat anti Maltose Phosphorylase antibody (200-101-230 lot 7850) was used to detect Maltose Phosphorylase under reducing (R) and non-reducing (NR) conditions.

Reduced samples of purified target proteins contained 4% BME and were boiled for 5 minutes. Samples of ~1µg of protein per lane were run by SDS-PAGE. Protein was transferred to nitrocellulose and probed with 1:3000 dilution of primary antibody (ON 4 C in MB-070). Detection shown was using Dylight 488 conjugated Donkey anti goat (605-741-125 lot 21094 1:10K in TBS/MB-070 1 hr RT).

Images were collected using the BioRad VersaDoc System

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.