

Datasheet for 200-106-212S**Bilirubin Oxidase Antibody Biotin Conjugated****Overview**

Description:	Anti-Bilirubin Oxidase (Myrothecium verrucaria) (GOAT) Antibody Biotin Conjugated - 200-106-212S
Item No.:	200-106-212S
Size:	25 µL
Applications:	WB
Reactivity:	Myrothecium verrucaria
Host Species:	Goat

Product Details

Background:	Anti-Bilirubin Oxidase antibody is an enzyme that increases the oxidation rate of bilirubin to biliverdin and other substrates. It is ideal for investigators interested in Metabolism or Signal Transduction research.
Synonyms:	goat anti-Bilirubin Oxidase biotin conjugation, goat anti-Bilirubin Oxidase Antibody biotin Conjugated
Host Species:	Goat
Conjugate:	Biotin
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	BLRO_MYRVE
Reactivity:	Myrothecium verrucaria
Immunogen Type:	Native Protein
Immunogen:	Bilirubin Oxidase [Myrothecium verrucaria]

Purity/Specificity: Anti-Bilirubin Oxidase is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Biotin, anti-Goat Serum as well as purified and partially purified Bilirubin Oxidase [Myrothecium verrucaria]. Cross reactivity against Bilirubin Oxidase from other sources is unknown.

Relevant Links:

- [UniProtKB - Q12737](#)

Application Details

Tested Applications: WB

Application Note: Anti-Bilirubin Oxidase has been tested by western blot and is suitable to be assayed against 1.0 ug of Bilirubin Oxidase in a standard capture ELISA using Peroxidase Conjugated Streptavidin #S000-03 and ABTS (2,2'-azino-bis-[3-ethylbenthiazoline-6-sulfonic acid]) code # ABTS-100 as a substrate for 30 minutes at room temperature. A working dilution of 1:1,000 to 1:3,000 of the reconstitution concentration is suggested for this product.

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

ELISA: 1:1,000 - 1:5,000

IP: 1:100

WB: 1:500 - 1:2,500

Formulation

Physical State: Liquid (sterile filtered)

Concentration: 1.0 mg/mL by UV absorbance at 280 nm

Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Preservative: 0.01% (w/v) Sodium Azide

Stabilizer: 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free


Shipping & Handling

Shipping Condition: Dry Ice

Storage Condition: Store vial at -20° C or below prior to opening. This vial contains a relatively low volume of reagent (25 µL). To minimize loss of volume dilute 1:10 by adding 225 µL of the buffer stated above directly to the vial. Recap, mix thoroughly and briefly centrifuge to collect the volume at the bottom of the vial. Use this intermediate dilution when calculating final dilutions as recommended below. Store the vial at -20°C or below after dilution. Avoid cycles of freezing and thawing.


Expiration: Expiration date is one (1) year from date of receipt.

Images



Western Blot

Western Blot of Goat anti-Bilirubin Oxidase (*Myrothecium verrucaria*) Antibody Biotin Conjugated. Lane 1: Bilirubin Oxidase (*Myrothecium verrucaria*). Load: 50 ng per lane. Primary antibody: Goat anti-Bilirubin Oxidase (*Myrothecium verrucaria*) Antibody Biotin Conjugated at 1:1,000 overnight at 4°C. Secondary antibody: HRP Streptavidin (p/n S000-03) 1:40,000 in MB-070 for 30 min at RT. Block: MB-070 for 30 minutes at RT. Predicted/Observed size: 64 kDa, 60 kDa for Bilirubin Oxidase. Other band(s): Bilirubin Oxidase splice variants and isoforms.



Western Blot

Rockland Goat anti Bilirubin Oxidase antibody (200-101-212 lot 6533) was used to detect Bilirubin Oxidase under reducing (R) conditions. A samples of ~1µg of purified protein contained 4% BME and was boiled for 5 minutes and run by SDS-PAGE on a 4-20% acrylamide gel. Protein was transferred to nitrocellulose and probed with 1:3000 dilution of primary antibody (ON 4 C in MB-070). Detection shown was using Dylight 488 conjugated Donkey anti goat (605-741-125 lot 21094 1:10K in TBS/MB-070 1 hr RT) and Dylight 649 conjugated Donkey anti goat (605-743-125 lot 20834 1:10K in TBS/MB-070) 1 hr RT. Images were collected using the BioRad VersaDoc System

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.