

Datasheet for 200-301-A37S**Wnt1 Antibody****Overview**

Description:	Anti-Wnt1 (MOUSE) Monoclonal Antibody - 200-301-A37S
Item No.:	200-301-A37S
Size:	25 µL
Applications:	WB
Reactivity:	Mouse
Host Species:	Mouse

Product Details**Background:**

The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. Wnt1 (Wingless-type MMTV integration site family member 1) is a member of the WNT gene family. It is highly conserved in evolution and the protein encoded by this gene is known to be 98% identical to mouse Wnt1 protein at the amino acid level. Studies in mouse indicate that the Wnt1 protein functions in the induction of the mesencephalon and cerebellum. This gene was originally considered as a candidate gene for Joubert syndrome, an autosomal recessive disorder with cerebellar hypoplasia as a leading feature. However, further studies suggested that the gene mutations might not have a significant role in Joubert syndrome. Wnt1 is secreted as an extracellular matrix protein.

Synonyms:	mouse anti-WNT1 antibody, INT1, Proto-oncogene protein Wnt-1
Host Species:	Mouse
Clonality:	Monoclonal
Clone ID:	13F9
Format:	IgG1

Target Details

Gene Name:	WNT1
Reactivity:	Mouse

Immunogen Type:	Conjugated Peptide
Immunogen:	Synthetic peptide corresponding to an internal region of human Wnt1 protein
Purity/Specificity:	Anti-Wnt1 antibody purified from tissue culture supernatant by Protein-A chromatography followed by extensive dialysis against the buffer stated above. This antibody reacts with human and mouse Wnt1 protein. A BLAST analysis was used to suggest cross-reactivity with Wnt1 from mouse, human, rat, bovine, dog, macaque, opossum and rat based on a 100% homology with the immunizing sequence. Partial cross-reactivity is expected against chicken Wnt1 based on a 91% sequence homology. Cross-reactivity with Wnt1 from other sources has not been determined.
Relevant Links:	• NCBI - 139743 • UniProtKB - P04628 • GeneID - 7471

Application Details

Tested Applications:	WB
Application Note:	This affinity purified antibody has been tested for use in western blotting. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 41 kDa in size corresponding to Wnt1 by western blotting in the appropriate cell lysate or extract.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:4,000 - 1:20,000
WB:	1:175 - 1:250

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1.0 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Sodium Azide
Stabilizer:	None


Shipping & Handling

Shipping Condition:	Dry Ice
----------------------------	---------

Storage Condition: Store vial at -20° C or below prior to opening. This vial contains a relatively low volume of reagent (25 µL). To minimize loss of volume dilute 1:10 by adding 225 µL of the buffer stated above directly to the vial. Recap, mix thoroughly and briefly centrifuge to collect the volume at the bottom of the vial. Use this intermediate dilution when calculating final dilutions as recommended below. Store the vial at -20°C or below after dilution. Avoid cycles of freezing and thawing.

Expiration: Expiration date is three (3) months from date of receipt.

Images



Western Blot

Western blot using Rockland's Protein A purified anti-Wnt1 monoclonal antibody shows detection of Wnt1 protein in mouse testis lysate. The results show specific binding corresponding to the ~41 kDa Wnt1 protein. Primary antibody was used at a 1:500 dilution. Personal communication, Stephen Brown, Brown University

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.