

Datasheet for 200-301-E58

Hsp40 (YDJ1) Antibody**Overview**

Description:	Anti-HSP40 (YDJ1) (MOUSE) Monoclonal Antibody - 200-301-E58
Item No.:	200-301-E58
Size:	100 µg
Applications:	IP, WB
Reactivity:	Yeast
Host Species:	Mouse

Product Details

Background:	Human Hsp40/DnaJ proteins comprise a large protein family, members of which feature the J domain (named after the bacterial DnaJ protein). The J-domain spans the first 75 N-terminal amino acids and is separated from the C-terminal by a glycine/phenylalanine-rich domain. There are two main types of Hsp40; type I DNAJ proteins including HDJ2 and yeast Ydj1; type II includes yeast Sis1 and human Hdj1. Whereas type I possesses a zinc finger domain which helps in the function of protein folding, type II does not. Members of the Hsp40/DnaJ family play diverse roles in many cellular processes, such as folding, translocation, degradation and assembly of multi-protein complexes. Hsp40 stimulates the ATPase activity of Hsp70 which in turn causes conformational changes of the unfolded proteins. The Hsp40-Hsp70-unfolded protein complex further binds to co-chaperones Hip, Hop and HSP90 which leads to protein folding, or components of protein degradation machinery CHIP and BAG-1.
Synonyms:	Mitochondrial protein import protein MAS5, Yeast dnaJ protein 1, YDJ1, MAS5, YNL064C, N2418, YNL2418C
Host Species:	Mouse
Clonality:	Monoclonal
Clone ID:	1G10.H8
Format:	IgG1

Target Details

Gene Name:	YDJ1
-------------------	------

Reactivity:	Yeast
Immunogen Type:	Recombinant Protein
Immunogen:	Hsp40 (YDJ1) Antibody was produced in mice by repeated immunizations raised against full length protein yeast Hsp40 (YDJ1).
Purity/Specificity:	Anti-Hsp40, YDJ1 Antibody was purified by Protein G chromatography. A BLAST analysis was used to suggest cross-reactivity with Hsp40 from yeast based on 100% homology with the immunizing sequence. Cross-reactivity with Hsp40 from other sources has not been determined. Heat Shock research.
Relevant Links:	<ul style="list-style-type: none">• NCBI - NP_014335.1• GeneID - 855661• UniProtKB - P25491

Application Details

Tested Applications:	IP, WB
Application Note:	Anti-Hsp40 Antibody has been tested by western blot and IP and is suitable for ELISA. Expect a band approximately ~40kDa protein corresponding to the molecular mass of Hsp40 on SDS PAGE immunoblots. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:200
IP:	User Optimized
WB:	1:2000

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1 mg/mL by UV absorbance at 280 nm
Preservative:	0.09% (w/v) Sodium Azide
Stabilizer:	50% (v/v) Glycerol


Shipping & Handling

Shipping Condition:	Dry Ice
----------------------------	---------

Storage Condition: Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Expiration: Expiration date is one (1) year from date of receipt.

Images



Western Blot

Western Blot of mouse anti-Hsp40 YDJ1 antibody. Lane 1: human cell line mix. Lane 2: none. Load: 10 µg per lane. Primary antibody: Hsp40 YDJ1 antibody at 1:1000 for overnight at 4°C. Secondary antibody: IRDye800™ mouse secondary antibody at 1:10,000 for 45 min at RT. Block: 5% BLOTTO overnight at 4°C. Predicted/Observed size: 44.7 kDa/~45 kDa for Hsp40 YDJ1. Other band(s): none.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.