

Datasheet for 200-301-E66**KLH Antibody****Overview**

Description:	Anti-KLH (MOUSE) Monoclonal Antibody - 200-301-E66
Item No.:	200-301-E66
Size:	100 µg
Reactivity:	KLH
Host Species:	Mouse

Product Details

Background:	KLH antibody is used to detect Keyhole Limpet Hemocyanin. KLH is a copper-containing oxygen carrier occurring freely dissolved in the hemolymph of many mollusks and arthropods. KLH is often used in molecular immunology as a carrier protein conjugated to low molecular weight molecules such as peptides, amino acids, nucleic acids, drugs or toxins to render them more immunogenic due to the size of the conjugate complex and the immunogenicity of KLH. KLH forms a large complex composed of ~50 kDa subunits.
Synonyms:	mouse anti-KLH Antibody, Keyhole Limpet Hemocyanin, KLH antibody
Host Species:	Mouse
Clonality:	Monoclonal
Clone ID:	3B1
Format:	IgM

Target Details

Reactivity:	KLH
Immunogen Type:	Native Protein
Immunogen:	Anti-KLH (MOUSE) Monoclonal Antibody was produced in mouse by repeated immunizations with KLH (Keyhole Limpet Hemocyanin) followed by hybridoma development.
Purity/Specificity:	Anti-KLH is a total IgM fraction prepared from ascites fluid by selective precipitation and tandem molecular sieve chromatography followed by extensive dialysis against the buffer stated above. KLH antibody is specific for KLH. Does not cross react with mouse heart or skin tissue sections. Cross-reactivity from other sources has not been determined.

Application Details

Application Note:	Anti-KLH can be used as an IgM isotype control in various assays or in ELISA assays to detect KLH or KLH conjugates. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:20,000 - 1:100,000
IHC:	1:1,000 - 1:5,000
WB:	1:2,000 - 1:10,000

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1.0 mg/mL by UV absorbance at 280 nm
Buffer:	0.1 M Tris Chloride, 0.5 M Sodium Chloride, pH 8.0
Preservative:	0.1% (w/v) Sodium Azide
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store KLH antibody at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.