

Datasheet for 200-301-E86

Alpha B Crystallin Antibody

Overview

Description:	Anti-Alpha B Crystallin (MOUSE) Monoclonal Antibody - 200-301-E86
Item No.:	200-301-E86
Size:	200 µg
Applications:	IF, IHC, WB
Reactivity:	Human, Bovine
Host Species:	Mouse

Product Details

Background:

The alpha-crystallins are major water-soluble lens structural proteins of the vertebrate eye that are related to the small heat shock protein family. The alpha-crystallins possess structural and functional similarities with Hsp25 and Hsp27. Mammalian lens crystallins are divided into alpha, beta and gamma families. Alpha and beta families are further divided into acidic and basic groups (Alpha-A and Alpha-B respectively). In the lens, alpha-crystallin primarily functions to maintain proper refractive index, however it can also function as a molecular chaperone that binds to the denatured proteins, keeping them in solution and thereby maintaining the translucency of the lens. When cellular stress occurs, alpha-crystallin enters its phosphorylated state and may serve a structural control function and play a role in protein maintenance. In addition to their interaction with proteins, alpha-crystallins also interact with native molecules such as membrane proteins, Golgi matrix protein, structural proteins, nuclear proteins and DNA. Two other functions are an autokinase activity and participation in the intracellular architecture, and it has also been proven that both alpha-A and B prevent apoptosis by inhibiting caspases. Specifically, alpha-B crystallin is found in many cells and organs outside the lens, and alpha B is overexpressed in several neurological disorders and in cell lines under stress conditions.

Synonyms:	AACRYA, CRYA2, CRYAB, CTPP2, HspB5, NY Ren 27 antigen antibody, alpha crystallin B antibody, Heat shock protein beta-5, Renal carcinoma antigen NY-REN-27, Rosenthal fiber component
Host Species:	Mouse
Clonality:	Monoclonal
Clone ID:	3A10.H4
Format:	IgG1

Target Details

Gene Name:	CRYAB
Reactivity:	Human, Bovine
Immunogen Type:	Native Protein
Immunogen:	Alpha B Crystallin Antibody was produced in mice by repeated immunizations with Native Alpha B Crystallin protein.
Purity/Specificity:	Anti-Alpha B Crystallin Antibody was purified by Protein G chromatography. A BLAST analysis was used to suggest cross-reactivity with Alpha B Crystallin from Human and bovine based on 100% homology with the immunizing sequence. Does not cross-react with α A-crystallin, β L-crystallin, BH-crystallin, γ -crystallin, Hsp25, Hsp27 or Hsp47 proteins. Cross-reactivity with Alpha B Crystallin from other sources has not been determined. Heat Shock research.
Relevant Links:	• NCBI - NP_001876.1 • GeneID - 1410 • UniProtKB - P02511

Application Details

Tested Applications:	IF, IHC, WB
Application Note:	Anti-Alpha B-crystallin Antibody has been tested in WB, IHC, and IF and is recommended for use in ELISA. Expect a band approximately ~20kDa (predicted mol. weight is ~21kDa) corresponding to α B-crystallin. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:200
WB:	1:2000


Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.09% (w/v) Sodium Azide
Stabilizer:	50% (v/v) Glycerol

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



Western Blot

Western Blot of mouse Anti-Alpha B Crystallin Antibody. Lane 1: Rat Brain. Lane 2: Rat Heart. Lane 3: Rat Kidney. Lane 4: Rat liver. Lane 5: Rat lung. Lane 6: Rat Pancreas. Lane 7: Rat skeletal muscle. Lane 8: Rat Spleen. Lane 9: Rat testes. Lane 10: Rat thymus. Lane 11: C2C12 cell line. Lane 12: NIH 3T3 ctrl. Lane 13: PAM212 cell line. Load: 20ug per lane. Primary antibody: AlphaB-Crystallin 1:1000 for overnight at 4°C. Secondary antibody: Goat anti-Mouse IgG HRP antibody at 1:40,000 for 45 min at RT. Block: 5% Blotto overnight at 4°C. Predicted/Observed size: 19.9 kDa, 20 kDa for Alpha B Crystallin.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.