

Datasheet for 200-301-F79**KChIP1 K⁺ channel Antibody****Overview**

Description:	Anti-KChIP1 K ⁺ channel (MOUSE) Monoclonal Antibody - 200-301-F79
Item No.:	200-301-F79
Size:	100 µg
Applications:	IHC
Reactivity:	Human, Mouse, Rat
Host Species:	Mouse

Product Details

Background:	There are four member of the KChIPs (Kv4 potassium channel interacting protein) family. They are all EF handcontaining proteins required for the traffic of channel forming Kv4 K ⁺ subunits to the plasma membrane. KChIP1 is expressed predominantly in the brain, with relative abundance in Purkinje cells of the cerebellum, the reticular neuclei, the medial habenular nuclei, the hippocampus and the striatum. It is expressed in a subpopulation of parvalbumin-positive neurons, suggesting its functional relationship with the GABAergic inhibitory neurons.
Synonyms:	Kcnp1, Kv channel-interacting protein 1, KChIP1, A-type potassium channel modulatory protein 1, Potassium channel-interacting protein 1
Host Species:	Mouse
Clonality:	Monoclonal
Clone ID:	S55-7
Format:	IgG1

Target Details

Gene Name:	Kcnp1
Reactivity:	Human, Mouse, Rat
Immunogen Type:	Recombinant Protein
Immunogen:	KChIP1 K ⁺ channel Antibody was produced in mice by repeated immunizations raised against a fusion protein of full length of rat KChIP1b.

Purity/Specificity: Anti-KChIP1 K⁺ channel Antibody was purified by Protein G chromatography. A BLAST analysis was used to suggest cross-reactivity with ChIP1 K⁺ channel from Mouse, human, and rat based on 100% homology with the immunizing sequence. No cross-reactivity against rat KChIPs2,3. Cross-reactivity with KChIP1 K⁺ channel from other sources has not been determined. Ion Channels research.

Relevant Links:

- [NCBI - NP_075218.2](#)
- [GeneID - 65023](#)
- [UniProtKB - Q8R426](#)

Application Details

Tested Applications:	IHC
Application Note:	Anti-KChIP1 K ⁺ channel Antibody is tested for use in IP and IHC. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
IHC:	0.1-1.0ug/mL


Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Stabilizer:	50% (v/v) Glycerol

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images

**Immunohistochemistry**

Immunohistochemistry of mouse anti-KChIP1 antibody.
Tissue: Frozen sections of mouse brain extract. Primary Antibody: KChIP1 antibody at 1:40 for 1h at RT. Secondary antibody: Peroxidase mouse secondary at 1:10,000 for 45 min at RT. Localization: cell membrane. Staining: KChIP1 as brown signal.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.