

Datasheet for 200-301-F87**Kir2.1 Antibody****Overview**

Description:	Anti-Kir2.1 (MOUSE) Monoclonal Antibody - 200-301-F87
Item No.:	200-301-F87
Size:	100 µg
Applications:	IF, IHC, IP, WB
Reactivity:	Human, Mouse, Rat
Host Species:	Mouse

Product Details

Background:	The Kir2.1 inward-rectifier potassium ion channel is encoded by the KCNJ2 gene. A defect in this gene is associated with Andersen-Tawil syndrome.
Synonyms:	HHBIRK1, HHIRK1, HIRK 1, IRK1, KCNJ2, LQT7, SQT3, potassium inwardly rectifying channel J2, Inward rectifier, potassium channel 2, Inward rectifier K(+) channel Kir2.1, IRK-1
Host Species:	Mouse
Clonality:	Monoclonal
Clone ID:	S112
Format:	IgG1

Target Details

Gene Name:	Kcnj2
Reactivity:	Human, Mouse, Rat
Immunogen Type:	Recombinant Protein
Immunogen:	Kir2.1 Antibody was produced in mice by repeated immunizations raised against a fusion protein at an internal region of mouse Kir2.1.
Purity/Specificity:	Anti-Kir2.1 Antibody was purified by Protein G chromatography. A BLAST analysis was used to suggest cross-reactivity with Kir2.1 from Human, Rat, and Mouse based on 100% homology with the immunizing sequence. No cross-reactivity against Kir2.2 or Kir2.3. Cross-reactivity with Kir2.1 from other sources has not been determined. Channels and Transporters research.

Relevant Links:	• NCBI - NP_032451 • GeneID - 16518 • UniProtKB - P35561
------------------------	--

Application Details

Tested Applications:	IF, IHC, IP, WB
Application Note:	Anti-Kir2.1 Antibody is tested for use in WB, IF, IP, and IHC-P. Expect a band approximately ~45kDa on specific lysates. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
IF:	1.0-10 µg/mL
IHC:	0.1-1.0 µg/mL
WB:	1.0 µg/mL


Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1.0 mg/mL by UV absorbance at 280 nm
Buffer:	1X PBS, pH 7.4
Preservative:	0.09% (w/v) Sodium Azide
Stabilizer:	50% (v/v) Glycerol

Shipping & Handling


Shipping Condition:	Dry Ice
Storage Condition:	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



Immunohistochemistry

Immunohistochemistry of mouse anti-Kir2.1 antibody.
 Tissue: Back sections of mouse. Primary Antibody: Kir2.1 antibody at 1 µg/mL for 1h at RT. Secondary antibody: Peroxidase mouse secondary at 1:10,000 for 45 min at RT. Localization: membrane. Staining: Kir2.1 as white signal.



Western Blot

Western Blot of mouse anti-Kir2.1 antibody. Lane 1: Cos Transient Cells. Primary antibody: Kir2.1 antibody at 1:1000 for overnight at 4°C. Secondary antibody: Goat anti-mouse IgG HRP secondary antibody at 1:10,000 for 45 min at RT. Block: 5% Blotto overnight 4°C. Predicted/Observed size: 48.3 kDa/55 kDa for Kir2.1 total. Other band(s): none.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.