

Datasheet for 200-301-F96**Kv1.6 Antibody****Overview**

Description:	Anti-Kv1.6 (MOUSE) Monoclonal Antibody - 200-301-F96
Item No.:	200-301-F96
Size:	100 µg
Applications:	IHC
Reactivity:	Human, Mouse, Rat
Host Species:	Mouse

Product Details

Background:	Voltage gated channels are tetrameters composed of four subunit consists of six transmembrane segments with cytoplasmic NH2 and COOH-termini. Members of the KV1-4 subfamilies generate functional K ⁺ channels in a homotetrameric configuration.
Synonyms:	HBK2, KCNA6, Potassium voltage-gated channel subfamily A member 6, MK1.6, Voltage-gated potassium channel subunit Kv1.6
Host Species:	Mouse
Clonality:	Monoclonal
Clone ID:	S19-36
Format:	IgG3

Target Details

Gene Name:	Kcna6
Reactivity:	Human, Mouse, Rat
Immunogen Type:	Conjugated Peptide
Immunogen:	Kv1.6 Antibody was produced in mice by repeated immunizations raised against a synthetic peptide at the cytoplasmic C-terminus region of rat Kv1.6.

Purity/Specificity: Anti-Kv1.6 Antibody was PEG purified. A BLAST analysis was used to suggest cross-reactivity with Kv1.6 from Mouse, Human, and rat based on 100% homology with the immunizing sequence. Cross-reactivity with Kv1.6 from other sources has not been determined. Ion Channels research.

Relevant Links:

- [UniProtKB - P17659](#)

Application Details

Tested Applications:	IHC
Application Note:	Anti-Kv1.6 Antibody is tested for use in IP and IHC. Expect a band approximately ~60kDa on specific lysates. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
IF:	1.0-10ug/mL
IHC:	0.1-1.0ug/mL
IP:	User Optimized
WB:	1ug/mL


Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.09% (w/v) Sodium Azide
Stabilizer:	50% (v/v) Glycerol

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



Immunohistochemistry

Immunohistochemistry of mouse anti-Kv1.6 antibody.
Tissue: Frozen sections of mouse brain extract. Primary Antibody: Kv1.6 antibody at 1 µg/mL for 1h at RT. Secondary antibody: Peroxidase mouse secondary at 1:10,000 for 45 min at RT. Localization: membrane. Staining: Kv1.6 as brown signal.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.