

Datasheet for 200-326-N63**CD11c Allophycocyanin Antibody****Overview**

Description:	Anti-CD11c (MOUSE) Allophycocyanin Conjugated Monoclonal Antibody - 200-326-N63
Item No.:	200-326-N63
Size:	500 µL
Applications:	FC
Reactivity:	Human
Host Species:	Mouse

Product Details

Background:	The 3.9 antibody is widely used as a marker for CD11c expression on dendritic cells (DC), often in parallel with markers for CD11b, for identification of developmental stages and mature subsets of this cell type. CD11c is prominently expressed on tissue macrophages, and is also detected on activated neutrophils, granulocytes, some types of activated T cells and intestinal intraepithelial lymphocytes (IEL).L killing through its interactions with fibrinogen, CD54, and iC3b.
Synonyms:	Integrin alpha-X, CD11 antigen-like family member C, Leu M5, Leukocyte adhesion glycoprotein p150,95 alpha chain, Leukocyte adhesion receptor p150 95, CD11c
Host Species:	Mouse
Conjugate:	Allophycocyanine (APC)
Clonality:	Monoclonal
Clone ID:	3.9
Format:	IgG1
F/P Ratio:	1-2

Target Details

Gene Name:	ITGAX
Reactivity:	Human
Immunogen:	Anti-CD11c Antibody (Monoclonal) was produced by repeated immunizations with CD11c antigen.

Purity/Specificity: Allophycocyanin conjugated CD11c Monoclonal Antibody was purified from tissue culture supernatant via affinity chromatography and is directed against human CD11c. Reactivity is observed against human CD11c, Baboon, Chimpanzee, Cynomolgus, and Rhesus. Cross reactivity with CD11c from other sources has not been tested. Anti-CD11c is conjugated with APC under optimal conditions and the solution is free of unconjugated APC.

Relevant Links:

- [UniProtKB - P20702](#)
- [NCBI - NP_000878.2](#)
- [GeneID - 3687](#)

Application Details

Tested Applications:	FC
Application Note:	Anti-CD11c is tested for Flow Cytometry and useful in Immunofluorescence and Immunoprecipitation. Researchers should determine optimal titers for applications that are not stated.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
FC:	5µL/test/1x10 ⁵ to 1x10 ⁸ cells
IF:	User Optimized
IP:	User Optimized

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	0.05 mg/mL Sufficient to run approximately 100 tests
Buffer:	0.01 M Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.09% (w/v) Sodium Azide
Stabilizer:	0.1% Gelatin


Shipping & Handling

Shipping Condition:	Wet Ice
Storage Condition:	Store vial at 4° C prior to opening. Dilute only prior to immediate use. This product is stable at 4° C as an undiluted liquid. Use subdued lighting during handling and incubation of cells prior to analysis. Store reagent in the dark. DO NOT FREEZE.

Expiration:


Expiration date is six (6) months from date of receipt.

Images



Flow Cytometry

Flow Cytometry of Mouse anti-CD11c Allophycocyanin Conjugated Monoclonal Antibody. Cells: human peripheral blood monocytes. Stimulation: none. Antibody: (Dotted Line) 0.25 μ g APC Mouse IgG1 isotype control; (BLUE) Allophycocyanin Anti-CD11c mouse antibody; (RED) Allophycocyanin Anti-CD11c mouse antibody using 5 μ l (0.25 μ g).



Flow Cytometry

Flow Cytometry of Mouse anti-CD11c Allophycocyanin Conjugated Monoclonal Antibody. Cells: human peripheral blood monocytes. Stimulation: none. Antibody: (Dotted Line) 0.25 μ g APC Mouse IgG1 isotype control; (BLUE) Allophycocyanin Anti-CD11c mouse antibody using 5 μ l (0.25 μ g).

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.