

Datasheet for 200-326-N76

CD154 (CD40L) Allophycocyanin Antibody**Overview**

Description:	Anti-CD154 (CD40L) (MOUSE) Allophycocyanin Conjugated Monoclonal Antibody - 200-326-N76
Item No.:	200-326-N76
Size:	500 µL
Applications:	FC
Reactivity:	Human
Host Species:	Mouse

Product Details

Background:	CD154 (CD40 ligand) is also known as CD40L, gp39, TRAP and T-BAM. CD40 ligand is a 32-39 kD type II transmembrane glycoprotein. It is a member of the TNF superfamily and is expressed on activated T cells. It has been reported to be important for B cell costimulation following binding of its receptor, CD40. Additionally, binding of CD40L to CD40 on B cells promotes the secretion of immunoglobulin and Ig isotype switching. CD40L is also involved in the regulation of cytokine production by T cells.
Synonyms:	CD40 ligand, CD40-L, T-cell antigen Gp39, TNF-related activation protein, TRAP, Tumor necrosis factor ligand superfamily member 5, CD40LG, CD40L, TNFSF5, TRAP
Host Species:	Mouse
Conjugate:	Allophycocyanine (APC)
Clonality:	Monoclonal
Clone ID:	Clone 24-31
Format:	IgG1
F/P Ratio:	1-2

Target Details

Gene Name:	CD40LG
Reactivity:	Human

Immunogen:	Anti-CD154 Antibody (Monoclonal) was produced by repeated immunizations with CD154 antigen.
Purity/Specificity:	Allophycocyanin conjugated CD154 Monoclonal Antibody was purified from tissue culture supernatant via affinity chromatography and is directed against human CD154. Reactivity is observed against human CD154. Cross reactivity is observed with Cynomolgus and Rhesus. Reactivity with CD154 from other sources has not been tested. Anti-CD154 is conjugated with APC under optimal conditions and the solution is free of unconjugated APC.
Relevant Links:	• UniProtKB - P29965 • NCBI - NP_000065.1 • GeneID - 959

Application Details

Tested Applications:	FC
Application Note:	Anti-CD154 is tested for Flow Cytometry and useful in Immunofluorescence. Researchers should determine optimal titers for applications that are not stated.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
FC:	5 uL/test/1x10 ⁵ to 1x10 ⁸ cells
IF:	User Optimized


Formulation

Physical State:	Liquid (sterile filtered)
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.09% (w/v) Sodium Azide
Stabilizer:	0.2% BSA (w/v)

Shipping & Handling

Shipping Condition:	Wet Ice
Storage Condition:	Store vial at 4° C prior to opening. Dilute only prior to immediate use. This product is stable at 4° C as an undiluted liquid. Use subdued lighting during handling and incubation of cells prior to analysis. Store reagent in the dark. DO NOT FREEZE.
Expiration:	Expiration date is six (6) months from date of receipt.

Images



Flow Cytometry

Flow Cytometry of Mouse anti-CD154 Allophycocyanin Conjugated Monoclonal Antibody. Cells: Human peripheral blood monocytes. Stimulation: TPA+ ionomycin (5 hrs). Antibody: (Unfilled Blue Line) Allophycocyanin Mouse IgG1 isotype control; (Blue Filled) Allophycocyanin Anti-CD154 mouse antibody using 5 μ L.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.