

Datasheet for 200-344-GQ1**ERK2 Internal Antibody DyLight™680****Overview**

Description:	Anti-ERK2 (MOUSE) Monoclonal Antibody DyLight™680 Conjugated - 200-344-GQ1
Item No.:	200-344-GQ1
Size:	100 µg
Applications:	WB
Reactivity:	Human, Mouse, Rat
Host Species:	Mouse

Product Details

Background:	ERK2 antibodies detect the ERK2 isoform. Mitogen activated protein kinase 1, also known as MAPK1, ERK, or ERK2, is an integral component of the MAP kinase cascade that regulates cell growth and differentiation. ERK1 and ERK2 are activated by MEK1 and MEK2 in the B-raf signaling pathway resulting in its translocation to the nucleus where it phosphorylates nuclear targets. Human ERK1 and ERK2 are 84% identical in sequence and share common functionality in cells. Anti-ERK2 antibody is ideal for investigators involved in Neuroscience, Cell Signaling and Cancer Research.
Synonyms:	mouse anti-ERK2 antibody DyLight™680 conjugation, DyLight™680 conjugated mouse anti-ERK2 antibody, MAPK1, ERK, ERK-2, ERK2, MAPK2, P42MAPK, PRKM1, PRKM2, DyLight™680, ERK-2 Antibody
Host Species:	Mouse
Conjugate:	DyLight™ 680
Clonality:	Monoclonal
Clone ID:	3H9.A1.A8
Format:	IgG
F/P Ratio:	1.54

Target Details

Gene Name:	MAPK1
-------------------	-------

Reactivity:	Human, Mouse, Rat
Immunogen Type:	Conjugated Peptide
Immunogen:	Anti-ERK2 Monoclonal Antibody was produced in mice by repeated immunizations with synthetic peptide corresponding to amino acid residues near an internal region of ERK2 conjugated to KLH.
Purity/Specificity:	This DyLight™680 conjugated protein A purified mouse monoclonal antibody reacts specifically with human ERK2. Anti-ERK2 is purified from tissue culture supernatant by protein A purification. Cross reactivity is expected to occur with human, mouse, and rat based on sequence identity of the peptide immunogen. This antibody does not react with the ERK1 isoform.
Relevant Links:	• UniProtKB - P28482

Application Details

Tested Applications:	WB
Application Note:	Anti-ERK 2 (MOUSE) DyLight™680 conjugated antibody is suitable for use in ELISA and Western Blotting. Specific conditions of reactivity should be optimized by the end user. Expect a band of approximately 41 kDa.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:40,000
WB:	1µg/ml


Formulation

Physical State:	Lyophilized
Concentration:	1mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Sodium Azide
Stabilizer:	10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free
Reconstitution Volume:	100 µL
Reconstitution Buffer:	Restore with deionized water (or equivalent)

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



Western Blot

Western Blot Results of Monoclonal Anti-ERK2 (MOUSE) Antibody DyLight™680 Conjugated.

Lane 1: Opal Prestained Molecular Weight Marker (p/n MB-210-0500).

Lane 2: human rERK1 Reduced [0.05µg].

Lane 3: ERK2 Reduced [0.05µg].

Primary Antibody: Anti-ERK2 (MOUSE) Antibody DyLight™680 Conjugated at 1.0µg/mL overnight at 2-8°C.

Secondary Antibody: Rabbit Anti-Mouse Antibody Peroxidase Conjugated (p/n 610-403-C46) at 1:40,000 for 30mins at RT.

Blocking buffer: Blocking Buffer for Fluorescent Western Blotting (p/n MB-070) for 1hr at RT.

Expected Molecular Weight: ~41kDa.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.