

**Datasheet for 200-345-215****GFP Antibody Dylight™ 800 Conjugated****Overview**

|                      |                                                                            |
|----------------------|----------------------------------------------------------------------------|
| <b>Description:</b>  | Anti-GFP (MOUSE) Monoclonal Antibody DyLight™ 800 Conjugated - 200-345-215 |
| <b>Item No.:</b>     | 200-345-215                                                                |
| <b>Size:</b>         | 100 µg                                                                     |
| <b>Applications:</b> | ELISA, WB                                                                  |
| <b>Reactivity:</b>   | GFP, eGFP, rGFP                                                            |
| <b>Host Species:</b> | Mouse                                                                      |

**Product Details**

|                      |                                                                                                                                                                                                                                                             |
|----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Background:</b>   | The emission spectra for this DyLight™ conjugate match the principle output wavelengths of most common fluorescence instrumentation.                                                                                                                        |
| <b>Synonyms:</b>     | mouse anti-GFP Antibody Dylight™ 800 Conjugation, Dylight™ 800 Conjugated mouse anti-GFP Antibody, GFP, Green Fluorescent Protein, GFP antibody, Green Fluorescent Protein antibody, EGFP, enhanced Green Fluorescent Protein, Aequorea victoria, Jellyfish |
| <b>Host Species:</b> | Mouse                                                                                                                                                                                                                                                       |
| <b>Conjugate:</b>    | DyLight™ 800                                                                                                                                                                                                                                                |
| <b>Clonality:</b>    | Monoclonal                                                                                                                                                                                                                                                  |
| <b>Clone ID:</b>     | 9F9.F9                                                                                                                                                                                                                                                      |
| <b>Format:</b>       | IgG1                                                                                                                                                                                                                                                        |
| <b>F/P Ratio:</b>    | 2.8                                                                                                                                                                                                                                                         |

**Target Details**

|                        |                                                                                                                                                                                              |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Reactivity:</b>     | GFP, eGFP, rGFP                                                                                                                                                                              |
| <b>Immunogen Type:</b> | Recombinant Protein                                                                                                                                                                          |
| <b>Immunogen:</b>      | BALB/c mice were immunized with a Green Fluorescent Protein (GFP) fusion protein corresponding to the full length amino acid sequence (246 aa) derived from the jellyfish Aequorea victoria. |

**Purity/Specificity:** GFP Dylight™ 800 Conjugated Antibody was prepared from tissue culture supernatant by Protein A affinity chromatography. Assay by Immunoelectrophoresis resulted in a single precipitin arc against anti-Mouse Serum. Reactivity is observed against recombinant Green Fluorescent Protein (000-001-215, recombinant GFP from Aequorea victoria) by Western blot. No reaction is seen against RFP.

**Relevant Links:**

- [UniProtKB - P42212](#)

## Application Details

**Tested Applications:** ELISA, WB

**Application Note:** Anti-GFP Dylight™ 800 Conjugated Antibody has been tested by ELISA and western blot and is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.

**Assay Dilutions:** All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

**FLISA:** >1:20,000

**IF:** >1:5,000

**WB:** 1:10,000 - 1:25,000

## Formulation

**Physical State:** Lyophilized

**Concentration:** 1.0 mg/mL by UV absorbance at 280 nm

**Buffer:** 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

**Preservative:** 0.01% (w/v) Sodium Azide

**Stabilizer:** 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free

**Reconstitution Volume:** 100 µL

**Reconstitution Buffer:** Restore with deionized water (or equivalent)


## Shipping & Handling

**Shipping Condition:** Ambient

**Storage Condition:** Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

**Expiration:** Expiration date is one (1) year from date of receipt.

## Images




### Western Blot

Two-color Western Blot using a DyLight™ 680 and DyLight™ 800 conjugates. Lane 1-5: 1μg, 0.5μg, 0.25μg, 0.125μg, and 0.063μg of Tubulin. Lane 6: Molecular weight. Lane 7-10: 6.2μg, 12.5μg, 2μg, 50μg of TNFα. Primary antibody: none. Secondary antibody: DyLight™ 680 and DyLight™ 800 mouse secondary antibody at 1:10,000. Block: MB-070 for 2 hrs at RT. Predicted size: 17kDa TNFα and 55kDa Tubulin. Other band(s): none.

### Diagram

Properties of DyLight™ Fluorescent Dyes. Fluorescence absorption and emission spectra of DyLight™ 800 in PBS, pH7.2. The emission spectra for this DyLight™ conjugate match the principle output wavelengths of most common fluorescence instrumentation.



## Western Blot

Properties of DyLight™ Conjugates.

| Emission      | Color                                                                             | DyLight™ Dye | Ex/Em (nm) | $\epsilon$ (M <sup>-1</sup> cm <sup>-1</sup> ) | Similar Dyes                       |
|---------------|-----------------------------------------------------------------------------------|--------------|------------|------------------------------------------------|------------------------------------|
| Blue          |  | 405          | 400/420    | 30,000                                         | Alexa™ 405, Cascade Blue           |
| Green         |  | 488          | 493/518    | 70,000                                         | Alexa™ 488, Cy2®, FITC             |
| Yellow        |  | 549          | 550/568    | 150,000                                        | Alexa™ 546, Alexa 555, Cy3®, TRITC |
| Red           |  | 649          | 646/674    | 250,000                                        | Alexa™ 647, Cy5®                   |
| Near Infrared |  | 680          | 682/715    | 140,000                                        | Alexa™ 680, Cy5.5®, IRDye™ 700     |
| Infrared      |  | 800          | 770/794    | 270,000                                        | IRDye™ 800                         |

## Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.