

Datasheet for 200-403-089-0100

## Dextranase Antibody Peroxidase Conjugated

### Overview

<b>Description:</b>	Anti-Dextranase (Penicillium species) (RABBIT) Antibody Peroxidase Conjugated - 200-403-089-0100
<b>Item No.:</b>	200-403-089-0100
<b>Size:</b>	100 µg
<b>Applications:</b>	ELISA, WB
<b>Reactivity:</b>	Penicillium
<b>Host Species:</b>	Rabbit

### Product Details

<b>Background:</b>	Dextranase belongs to the glycosyl hydrolase 49 family. It's catalytic reaction is endohydrolysis of (1->6)-alpha-D-glucosidic linkages in dextran.
<b>Synonyms:</b>	rabbit anti-Dextranase Antibody Peroxidase Conjugation, HRP Conjugated rabbit anti-Dextranase Antibody, dextranase, Alpha-1,6-glucan-6-glucanohydrolase, DEX
<b>Host Species:</b>	Rabbit
<b>Conjugate:</b>	Peroxidase (HRP)
<b>Clonality:</b>	Polyclonal
<b>Format:</b>	IgG

### Target Details

<b>Gene Name:</b>	DEX
<b>Reactivity:</b>	Penicillium
<b>Immunogen Type:</b>	Native Protein
<b>Immunogen:</b>	Dextranase [Penicillium species]

**Purity/Specificity:** Dextranase is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Rabbit Serum as well as purified and partially purified Dextranase [Penicillium species]. Cross reactivity against Dextranase from other sources is unknown.

**Relevant Links:**

- [UniProtKB - P48845](#)
- [NCBI - AAB47720.1](#)

## Application Details

**Tested Applications:** ELISA, WB

**Application Note:** Anti-Dextranase Peroxidase has been tested by ELISA and western blot and assayed against 1.0 ug of Dextranase [Penicillium species] in a standard capture ELISA using ABTS (2,2'-azino-bis-[3-ethylbenzothiazoline-6-sulfonic acid]) code # ABTS-100 as a substrate for 30 minutes at room temperature. A working dilution of 1:1,000 to 1:3,000 of the reconstitution concentration is suggested for this product.

**Assay Dilutions:** All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

**ELISA:** 1:2,000 - 1:10,000

**IP:** User Optimized

**WB:** 1:500 - 1:5,000

## Formulation

**Physical State:** Lyophilized

**Concentration:** 1.0 mg/mL by UV absorbance at 280 nm

**Buffer:** 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

**Preservative:** 0.01% (w/v) Gentamicin Sulfate. Do NOT add Sodium Azide!

**Stabilizer:** 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free

**Reconstitution Volume:** 100 µL

**Reconstitution Buffer:** Restore with deionized water (or equivalent)

## Shipping & Handling


**Shipping Condition:** Ambient

---

<b>Storage Condition:</b>	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
<b>Expiration:</b>	Expiration date is one (1) year from date of receipt.


---

## Images



### Western Blot

Western Blot of Rabbit Anti-Dextranase Antibody Peroxidase Conjugated. Lane 1: Opal Prestained Molecular Weight Marker (p/n MB-210-0500). Lane 2: Dextranase. Primary Antibody: Rabbit Anti-Dextranase at 1:1000 overnight at 2-8°C. Secondary Antibody: none. Block: MB-070. Expect MW: ~66kDa.



### ELISA

ELISA Results of Rabbit Anti-Dextranase Antibody Peroxidase Conjugate tested against purified Dextranase. Each well was coated in duplicate with 1.0 µg of Dextranase. The starting dilution of antibody was 5µg/ml and the X-axis represents the Log10 of a 3-fold dilution. This titration is a 4-parameter curve fit where the IC50 is defined as the titer of the antibody. Assay performed using HRP Conjugate Stabilizer (p/n MB-076) and TMB substrate (p/n TMBE-1000).

## Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.