

Datasheet for 200-403-140S**Adenosine Deaminase Antibody Peroxidase Conjugated****Overview**

Description:	Anti-Adenosine Deaminase (Calf Spleen) (RABBIT) Antibody Peroxidase Conjugated - 200-403-140S
Item No.:	200-403-140S
Size:	25 µL
Applications:	ELISA, WB
Reactivity:	Bovine Calf Spleen
Host Species:	Rabbit

Product Details

Background:	Adenosine deaminase (adenosine aminohydrolase, or ADA) is an enzyme involved in purine metabolism. It is needed for the breakdown of adenosine from food and for the turnover of nucleic acids in tissues. ADA irreversibly deaminates adenosine, converting it to the related nucleoside inosine by the substitution of the amino group for a hydroxyl group. ADA irreversibly deaminates adenosine, converting it to the related nucleoside inosine by the substitution of the amino group for a hydroxyl group. Inosine can then be dephosphorylated by another enzyme called purine nucleoside phosphorylase (PNP), converting it to hypoxanthine.
Synonyms:	rabbit anti-Adenosine Deaminase Antibody Peroxidase Conjugation, HRP conjugated rabbit anti-Adenosine Deaminase Antibody, Adenosine deaminase, Adenosine aminohydrolase
Host Species:	Rabbit
Conjugate:	Peroxidase (HRP)
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	ADA
Reactivity:	Bovine Calf Spleen
Immunogen Type:	Native Protein

Immunogen:	Adenosine Deaminase [Calf Spleen] ADA
Purity/Specificity:	Adenosine deaminase is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Rabbit Serum as well as purified and partially purified Adenosine Deaminase [Calf Spleen]. Cross reactivity against Adenosine Deaminase from other sources is unknown.
Relevant Links:	• UniProtKB - P56658 • NCBI - NP_776312.1 • GenelD - 280712

Application Details

Tested Applications:	ELISA, WB
Application Note:	Anti-Adenosine deaminase Peroxidase antibody has been tested by ELISA and western blot and is assayed against 1.0 ug of Adenosine Deaminase [Calf Spleen] in a standard capture ELISA using ABTS (2,2'-azino-bis-[3-ethylbenzothiazoline-6-sulfonic acid]) code # ABTS-100 as a substrate for 30 minutes at room temperature. A working dilution of 1:20,000 to 1:100,000 of the reconstitution concentration is suggested for this product.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:1,000 - 1:5,000
IP:	1:100
WB:	1:500 - 1:5,000

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1.0 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Gentamicin Sulfate. Do NOT add Sodium Azide!
Stabilizer:	10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free


Shipping & Handling

Shipping Condition: Dry Ice

Storage Condition: Store vial at -20° C or below prior to opening. This vial contains a relatively low volume of reagent (25 µL). To minimize loss of volume dilute 1:10 by adding 225 µL of the buffer stated above directly to the vial. Recap, mix thoroughly and briefly centrifuge to collect the volume at the bottom of the vial. Use this intermediate dilution when calculating final dilutions as recommended below. Store the vial at -20°C or below after dilution. Avoid cycles of freezing and thawing.


Expiration: Expiration date is one (1) year from date of receipt.

Images



Western Blot

Rockland biotin conjugated anti adenosine deaminase (200-406-140) was used to detect adenosine deaminase in mouse pancreas lysate (Left, Lane 1, 30 ul) under reducing conditions. The antibody was also used to detect purified Adenosine Deaminase (right, Lane 2), and endogenous Adenosine Deaminase in whole cell lysate from Jurkat and Raji cells (1:1 mixture, lane 3,) as well as Mouse Pancreas and Liver (1:1 mixture, lane 4). Lysates were run on 4-20% gel 140V under reducing conditions, transferred for 30 minutes at 100 V and blocked with 3% Fish Gel (left) or 3% BSA. Blot was incubated with 200-406-140 lot 5628 (1:5K in TBS, ON 4°C), washed 3X in TBS and incubated for 30 minutes with Dylight 488 conjugated Streptavidin (S000-41 lot 20833 1:5K in MB-070). Blot was imaged with the Biorad VersaDoc imaging system.



Western Blot

Western Blot of Rabbit anti-Adenosine Deaminase Antibody Peroxidase Conjugated. Lane 1: Adenosine Deaminase (Calf Spleen). Load: 50 ng per lane. Primary antibody: Rabbit anti-Adenosine Deaminase Antibody Peroxidase Conjugated at 1:1,000 overnight at 4°C. Secondary antibody: n/a. Block: MB-070 for 30 min at RT. Predicted/Observed size: 41 kDa, 45 kDa for Adenosine Deaminase (Calf Spleen).

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.