

Datasheet for 200-406-142S**ALKALINE PHOSPHATASE Antibody Biotin Conjugated****Overview**

Description:	Anti-Alkaline Phosphatase (Human Intestine) (RABBIT) Antibody Biotin Conjugated - 200-406-142S
Item No.:	200-406-142S
Size:	25 µL
Applications:	Dot Blot, ELISA
Reactivity:	Human
Host Species:	Rabbit

Product Details

Background:	At a high pH, Anti-Alkaline Phosphatase antibodies strip phosphate groups from DNA, RNA, and proteins. ALP is most abundant in placental and intestinal tissue specific such as the liver, placenta, and bone. Located in the cell membrane, Anti-Alkaline Phosphatase antibodies are ideal for researchers involved in Cancer, Enzyme, Tags and Cell Marker, and Cell Cycle research.
Synonyms:	rabbit anti-Alkaline Phosphatase Antibody biotin Conjugation, biotin conjugated rabbit anti-Alkaline phosphatase intestinal antibody, Alkaline phosphomonoesterase antibody, ALPI antibody, Glycerophosphatase antibody, IAP antibody, Intestinal alkaline phosphatase antibody, Kasahara isozyme antibody
Host Species:	Rabbit
Conjugate:	Biotin
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	ALP1
Reactivity:	Human
Immunogen Type:	Native Protein
Immunogen:	Alkaline Phosphatase [Human Intestine]

Purity/Specificity: Anti-Alkaline Phosphatase is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Biotin, anti-Rabbit Serum as well as purified and partially purified Alkaline Phosphatase [Human Intestine]. Cross reactivity against Alkaline Phosphatase from other sources may occur but has not been specifically determined.

Relevant Links:

- [UniProtKB - P09923](#)
- [NCBI - NP_001622.2](#)
- [GeneID - 248](#)

Application Details

Tested Applications: Dot Blot, ELISA

Application Note: Anti-Alkaline Phosphatase Biotin Antibody has been tested by ELISA and dot blot. This product is assayed against 1.0 ug of Alkaline Phosphatase in a standard capture ELISA using Peroxidase Conjugated Streptavidin #S000-03 and ABTS (2,2'-azino-bis-[3-ethylbenthiazoline-6-sulfonic acid]) code # ABTS-100 as a substrate for 30 minutes at room temperature. A working dilution of 1:4,000 to 1:20,000 of the reconstitution concentration is suggested for this product.

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

ELISA: User Optimized

WB: User Optimized

Formulation

Physical State: Liquid (sterile filtered)

Concentration: 1 mg/mL by UV absorbance at 280 nm

Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Preservative: 0.01% (w/v) Sodium Azide

Stabilizer: 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free


Shipping & Handling

Shipping Condition: Dry Ice

Storage Condition: Store vial at -20° C or below prior to opening. This vial contains a relatively low volume of reagent (25 µL). To minimize loss of volume dilute 1:10 by adding 225 µL of the buffer stated above directly to the vial. Recap, mix thoroughly and briefly centrifuge to collect the volume at the bottom of the vial. Use this intermediate dilution when calculating final dilutions as recommended below. Store the vial at -20°C or below after dilution. Avoid cycles of freezing and thawing.

Expiration: Expiration date is one (1) year from date of receipt.

Images



Dot Blot

Dot Blot for Anti-ALKALINE PHOSPHATASE (Human Intestine) (RABBIT) Antibody Biotin Conjugated. Lane 1: 100 ng Alkaline Phosphatase. Lanes 2-5: serial dilution 3 fold of Antigen. Primary Antibody: Anti-Alkaline Phosphatase Antibody Biotin Conjugated at 10 µg/ml at RT for 1 hour. Secondary Antibody: Streptavidin Peroxidase Conjugated (p/n S000-03) at a 1:40,000 dilution at RT for 30 minutes. Block: MB-070 at RT for 30 minutes.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.