

Datasheet for 200-4340**Acid Phosphatase Antibody Peroxidase Conjugated****Overview**

Description:	Anti-Acid Phosphatase (Potato) (RABBIT) Antibody Peroxidase Conjugated (BULK ORDER) - 200-4340
Item No.:	200-4340
Size:	20 mg
Applications:	Dot Blot, WB
Reactivity:	Acid Phosphatase (Potato)
Host Species:	Rabbit

Product Details

Background:	The catalytic activity of Acid Phosphatase is phosphate monoester + H ₂ O = an alcohol + phosphate.
Synonyms:	rabbit anti-Acid Phosphatase Antibody Peroxidase Conjugation, HRP conjugated rabbit anti-Acid Phosphatase Antibody, Acid phosphatase 1 soluble antibody, Acid phosphatase of erythrocyte antibody, Adipocyte acid phosphatase antibody, Cytoplasmic phosphotyrosyl protein phosphatase antibody, HAAP antibody
Host Species:	Rabbit
Conjugate:	Peroxidase (HRP)
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	PAP3
Reactivity:	Acid Phosphatase (Potato)
Immunogen Type:	Native Protein
Immunogen:	Acid Phosphatase [Potato]

Purity/Specificity:	This product is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-peroxidase, anti-Rabbit Serum as well as purified and partially purified Acid Phosphatase [Potato]. Cross reactivity against Acid Phosphatase from other tissues and species may occur but have not been specifically determined.
----------------------------	---

Relevant Links:	• UniProtKB - Q6J5M8 • NCBI - AAT37528.1
------------------------	--

Application Details

Tested Applications:	Dot Blot, WB
Application Note:	Anti-Acid Phosphatase Peroxidase Conjugated Antibody has been tested by dot blot and western blot and is suitable for use in ELISA and immunohistochemistry. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:4,000 - 1:20,000
WB:	1:500 - 1:2,000


Formulation

Physical State:	Lyophilized
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Gentamicin Sulfate. Do NOT add Sodium Azide!
Reconstitution Volume:	2.0 mL
Reconstitution Buffer:	Restore with deionized water (or equivalent)

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



Western Blot

Western Blot of Acid Phosphatase (Potato) Antibody. Lane 1: Acid Phosphatase (Potato). Load: 50 ng per lane. Primary antibody: Acid Phosphatase Antibody at 1:1,000 overnight at 4°C. Secondary antibody: HRP rabbit secondary antibody at 1:40,000 for 30 min at RT. Block: MB-070 for 30 min at RT. Predicted/Observed size: 49 kDa/42-45kDa for Acid Phosphatase (Potato) Protein. Other band(s): Acid Phosphatase monomer fragmentation.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.