

Datasheet for 200-4368-0100**Elastase Antibody Peroxidase Conjugated****Overview**

Description:	Anti-ELASTASE (RABBIT) Antibody Peroxidase Conjugated - 200-4368-0100
Item No.:	200-4368-0100
Size:	100 µg
Applications:	ELISA, WB
Reactivity:	Pig
Host Species:	Rabbit

Product Details

Background:	Elastase is an enzyme from the class of proteases (peptidases) that break down proteins. Elastase breaks down elastin, an elastic fiber that, together with collagen, determines the mechanical properties of connective tissue. Elastase also has the important immunological role of breaking down Shigella virulence factors. Elastase is responsible for pulmonary emphysema, cyclic hematopoiesis, and blistering in bullous pemphigoid.
Synonyms:	rabbit anti-Elastase Antibody Peroxidase Conjugation, HRP Conjugated rabbit anti-Elastase Antibody, Fecal pancreatic elastase 1 antibody, HLE antibody, Medullasin antibody, Pancreatic elastase IIA antibody, PE 1 antibody, PE1 antibody, Protease E antibody
Host Species:	Rabbit
Conjugate:	Peroxidase (HRP)
Clonality:	Polyclonal
Format:	IgG
F/P Ratio:	0.455

Target Details

Gene Name:	CELA1
Reactivity:	Pig
Immunogen Type:	Native Protein

Immunogen:	Elastase [Porcine Pancreas]
Purity/Specificity:	Elastase is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Rabbit Serum as well as purified and partially purified Elastase [Porcine Pancreas]. Cross reactivity against Elastase from other tissues and species may occur but have not been specifically determined.
Relevant Links:	• UniProtKB - P00772 • NCBI - NP_998988.1 • GenelD - 396766

Application Details

Tested Applications:	ELISA, WB
Application Note:	Anti-Elastase Peroxidase conjugated antibody has been tested by ELISA and western blot. This product is assayed against 1.0 ug of Elastase [Porcine Pancreas] in a standard capture ELISA using ABTS (2,2'-azino-bis-[3-ethylbenzthiazoline-6-sulfonic acid]) code # ABTS-100 as a substrate for 30 minutes at room temperature. A working dilution of 1:200 to 1:1,000 of the reconstitution concentration is suggested for this product.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:5,000 - 1:20,000
IHC:	User Optimized
WB:	1:500 - 1:2,000

Formulation

Physical State:	Lyophilized
Concentration:	1.0 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Gentamicin Sulfate. Do NOT add Sodium Azide!
Stabilizer:	10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free
Reconstitution Volume:	100 µL
Reconstitution Buffer:	Restore with deionized water (or equivalent)

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.