

Datasheet for 200-4675S**Penicillinase Antibody Biotin Conjugated****Overview**

Description:	Anti-Penicillinase (Enterobacter cloacae) (RABBIT) Antibody Biotin Conjugated - 200-4675S
Item No.:	200-4675S
Size:	25 µL
Applications:	Dot Blot
Reactivity:	Enterobacter cloacae
Host Species:	Rabbit

Product Details

Background:	Anti-Penicillinase antibody recognizes the protein penicillinase, a beta lactamase specific to penicillins. Penicillinase hydrolyzes the beta-lactam ring of penicillins destroying their antibacterial potency. Thus, penicillinase is often found in antibiotic-resistant bacteria.
Synonyms:	rabbit anti-Penicillinase Antibody biotin Conjugation, biotin Conjugated rabbit anti-Penicillinase Antibody, Aer1 antibody, Beta lactamase 1 antibody, Beta lactamase AER 1 antibody, OXA9 antibody, Oxacillinase antibody, PenPC antibody, Cephalosporinase
Host Species:	Rabbit
Conjugate:	Biotin
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	ampC
Reactivity:	Enterobacter cloacae
Immunogen Type:	Native Protein
Immunogen:	Penicillinase [Enterobacter cloacae]

Purity/Specificity: Anti-Penicillinase antibody is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Biotin, anti-Rabbit Serum as well as purified and partially purified Penicillinase [*Enterobacter cloacae*]. Cross reactivity against Penicillinase from other sources is unknown.

Relevant Links:

- [UniProtKB - P05364](#)
- [NCBI - AA024600.1](#)
- [GeneID - 9125732](#)

Application Details

Tested Applications: Dot Blot

Application Note: Anti-Penicillinase Biotin Conjugated antibody has been tested by dot blot and is suitable to be assayed against 1.0 ug of Penicillinase in a standard capture ELISA using Peroxidase Conjugated Streptavidin #S000-03 and ABTS (2,2'-azino-bis-[3-ethylbenthiazoline-6-sulfonic acid]) code # ABTS-100 as a substrate for 30 minutes at room temperature. A working dilution of 1:4,000 to 1:20,000 of the reconstitution concentration is suggested for this product.

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

IP: User Optimized

WB: User Optimized

Formulation

Physical State: Liquid (sterile filtered)

Concentration: 1.0 mg/mL by UV absorbance at 280 nm

Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Preservative: 0.01% (w/v) Sodium Azide

Stabilizer: 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free

Shipping & Handling

Shipping Condition: Dry Ice

Storage Condition: Store vial at -20° C or below prior to opening. This vial contains a relatively low volume of reagent (25 µL). To minimize loss of volume dilute 1:10 by adding 225 µL of the buffer stated above directly to the vial. Recap, mix thoroughly and briefly centrifuge to collect the volume at the bottom of the vial. Use this intermediate dilution when calculating final dilutions as recommended below. Store the vial at -20°C or below after dilution. Avoid cycles of freezing and thawing.

Expiration: Expiration date is one (1) year from date of receipt.

Images



Dot Blot

Dot Blot results of Rabbit Anti-Penicillinase Antibody Biotin Conjugated. Dots are Penicillinase at (1) 100ng, (2) 33.3ng, (3) 11.1ng, (4) 3.70ng, (5) 1.23ng. Blocking: MB-070 for 30 min at RT. Primary Antibody: Rabbit Anti-Penicillinase Biotin at 1µg/mL for 1hr at RT. Secondary Antibody: Streptavidin-HRP (p/n S000-03) at 1:40,000 for 30min at RT. Imaged with BioRad ChemiDoc, Chemi filter.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.