

Datasheet for 200-501-I26**TLR4 Antibody****Overview**

Description:	Anti-TLR4 (RAT) Monoclonal Antibody - 200-501-I26
Item No.:	200-501-I26
Size:	100 µg
Applications:	FC, IF, IP
Reactivity:	Mouse
Host Species:	Rat

Product Details

Background: Anti-TLR4 CD284 Antibody detects human TLR4. The Toll-like receptor (TLR) family in mammal comprises a family of transmembrane proteins characterized by multiple copies of leucine rich repeats in the extracellular domain and IL-1 receptor motif in the cytoplasmic domain. Like its counterparts in Drosophila, TLRs signal through adaptor molecules. The TLR family is a phylogenetically conserved mediator of innate immunity that is essential for microbial recognition. Ten human homologs of TLRs (TLR1-10) have been described. Among this family of receptors, TLR2 and TLR4 have been most studied. These studies have suggested that TLR2 and TLR4 may serve as potential main mediators of LPS signaling. The TLR4 cDNA codes for a protein consisting of 799 amino acids with approximate molecular weight of 88 kDa. Anti-TLR4 CD284 Antibody is ideal for investigators involved in cytokines and growth factor research.

Synonyms:	hToll, Toll-like receptor 4, CD284
Host Species:	Rat
Clonality:	Monoclonal
Clone ID:	MTS510
Format:	IgG2a

Target Details

Gene Name:	Tlr4
Reactivity:	Mouse

Immunogen Type:	Other
Immunogen:	TLR4 Antibody was produced in rats prepared by repeated immunizations with the Ba/F3 cell line expressing mouse protein TLR4.
Purity/Specificity:	Anti-TLR4 Antibody was purified by Protein G chromatography. A BLAST analysis was used to suggest cross-reactivity with Anti-TLR4 from mouse based on 100% homology with the immunizing sequence. Cross-reactivity with Anti-TLR4 from other sources has not been determined. Toll Like Receptors and signaling research.
Relevant Links:	<ul style="list-style-type: none">• UniProtKB - Q9QUK6• NCBI - NP_067272.1• GeneID - 21898

Application Details

Tested Applications:	FC, IF, IP
Application Note:	Anti-TLR4 Antibody has been tested in ICC/IF, Flow (CS), and IP. Expect a band approximately 95kDa on specific lysates. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
FC:	2-5 µg/10 ⁶ cells
IP:	User Optimized
WB:	1:500-1:1000

Formulation


Physical State:	Liquid
Concentration:	0.5 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.05% (w/v) Sodium Azide
Stabilizer:	0.05% BSA

Shipping & Handling

Shipping Condition:	Dry Ice
----------------------------	---------

Storage Condition:	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



Flow Cytometry

Flow Cytometry of Rat Anti-TLR4 CD284 antibody. Cells: mouse peritoneal. Stimulation: none. Primary Antibody: TLR4 CD284 antibody at 1 µg/1x10⁶ cells (red) and isotype control (green).

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.