

Datasheet for 200-601-144-0100**Aspartate Transaminase Antibody****Overview**

Description:	Anti-Aspartate Transaminase (AST) (SHEEP) Antibody - 200-601-144-0100
Item No.:	200-601-144-0100
Size:	100 µg
Reactivity:	Pig
Host Species:	Sheep

Product Details

Background:	Anti-Aspartate Transaminase is detects Aspartate Transaminase. Aspartate Transaminase is a pyridoxal phosphate (PLP)-dependent transaminase enzyme. AST catalyzes the reversible transfer of an α -amino group between aspartate and glutamate and, as such, is an important enzyme in amino acid metabolism. Anti-Aspartate Transaminase antibody is ideal for investigators involved in Cell Signaling, cellular biology and Signal Transduction research.
Synonyms:	sheep anti-Aspartate Transaminase Antibody, Aspartate aminotransferase 1 antibody, Aspartate aminotransferase cytoplasmic antibody, Aspartate aminotransferase cytosolic antibody, ec 2.6.1.1 antibody, GIG 18 antibody, GIG18 antibody, Glutamate oxaloacetate transaminase 1 antibody, AST
Host Species:	Sheep
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	GOT1
Reactivity:	Pig
Immunogen Type:	Native Protein
Immunogen:	Anti-Aspartate Transaminase antibody was produced by repeated immunizations with Aspartate Transaminase from pig heart.

Purity/Specificity: Anti-Aspartate Transaminase is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Sheep Serum as well as purified and partially purified Aspartate Transaminase (AST) [Pig Heart]. Cross reactivity against Aspartate Transaminase (AST) from other tissues and species may occur but have not been specifically determined.

Relevant Links:

- [UniProtKB - P00503](#)
- [NCBI - NP_999092.1](#)
- [GeneID - 396967](#)

Application Details

Application Note: Anti-Aspartate Transaminase is suitable for western blotting and for ELISA. Researchers should determine optimal titers for applications that are not stated below.

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

ELISA: 1:5,000 - 1:25,000

WB: 1:500 - 1:3,000

Formulation

Physical State: Lyophilized

Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Preservative: 0.01% (w/v) Sodium Azide

Stabilizer: None

Reconstitution Volume: 100 µL

Reconstitution Buffer: Restore with deionized water (or equivalent)

Shipping & Handling

Shipping Condition: Ambient

Storage Condition: Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Expiration: Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.