

Datasheet for 200-901-EF4**SDHAF2 Antibody****Overview**

Description:	Anti-SDHAF2 (CHICKEN) Antibody - 200-901-EF4
Item No.:	200-901-EF4
Size:	100 µg
Applications:	ELISA, WB
Reactivity:	Human, Mouse, Rat
Host Species:	Chicken

Product Details

Background:	SDHAF2 (Succinate dehydrogenase complex assembly factor 2) encodes a mitochondrial protein needed for FAD cofactor attach to the SDH enzyme which plays a critical role in mitochondria. SDHAF2 is a tumor suppressor, but recent evidence suggests that somatic mutations of the SDHAF2 are unlikely to contribute to parathyroid tumor development in sporadic primary hyperparathyroidism.
Synonyms:	SDHAF2 Antibody, PGL2, SDH5, C11orf79, PGL2, Succinate dehydrogenase assembly factor 2, mitochondrial, Succinate dehydrogenase subunit 5, SDH assembly factor 2
Host Species:	Chicken
Clonality:	Polyclonal
Format:	IgY

Target Details

Gene Name:	SDHAF2
Reactivity:	Human, Mouse, Rat
Immunogen Type:	Conjugated Peptide
Immunogen:	Anti-SDHAF2 antibody was prepared from chicken egg fractions produced by repeated immunizations with a 16 amino acid synthetic peptide near the N-terminus of human SDHAF2.

Purity/Specificity:	Anti-SDHAF2 Antibody is an IgY fraction that has been affinity purified by immunoaffinity purification. SDHAF2 antibody is predicted to not cross-react with other SDHAF protein family members.
Relevant Links:	• UniProtKB - Q9NX18 • GeneID - 54949 • NCBI - NP_060311

Application Details

Tested Applications:	ELISA, WB
Application Note:	Anti-SDHAF2 Antibody has been tested for use in ELISA and Western Blotting. Specific conditions for reactivity should be optimized by the end user. Expect a band at approximately 20 kDa in Western Blots of specific cell lysates and tissues.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:10,000
WB:	1 µg/mL


Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1 mg/mL by UV absorbance at 280 nm
Buffer:	0.01 M Sodium Phosphate, 0.25 M Sodium Chloride, pH 7.2
Preservative:	0.02% (w/v) Sodium Azide
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



Western Blot

Western blot of SDHAF2.

Load: rat liver tissue lysate.

Primary Antibody: SDHAF2 antibody at 1 μ g/mL in (A) the absence and (B) the presence of blocking peptide.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.