

Datasheet for 200-B02-N92**TCRbeta Fluorescein Antibody****Overview**

Description:	Anti-TCRbeta (Armenian Hamster) Fluorescein Conjugated Monoclonal Antibody - 200-B02-N92
Item No.:	200-B02-N92
Size:	500 µg
Applications:	FC
Reactivity:	Mouse
Host Species:	Hamster Armenian

Product Details

Background: The H57-597 antibody is specific for the beta chain of the mouse T cell Receptor (TCR). This cell surface protein combines with a second protein chain (alpha chain) to form the alpha-beta TCR that is expressed by NK1.1+ thymocytes, NKT cells, and the majority of peripheral T cells. A small number of T cells may express an alternative heteromer of gamma and delta protein chains, known as the gamma-delta TCR. These receptors participate in a complex with CD3, and with the co-receptors CD4 or CD8, to recognize and respond to antigens bound to MHC molecules on antigen-presenting cells. Such interactions promote T cell receptor signaling (T cell activation) and can result in a number of cellular responses including proliferation, differentiation, production of cytokines or activation-induced cell death. The H57-597 antibody is used as a phenotypic marker for T cells expressing the alpha-beta TCR. It is also widely used to cross-link surface TCR and thereby mimic TCR-mediated cell activation or induction of apoptosis. The antibody does not cross-react with cells expressing the gamma-delta TCR.

Synonyms:	Tib, TCRbeta, T cell receptor beta chain, Tcrb
Host Species:	Hamster Armenian
Conjugate:	Fluorescein (FITC)
Clonality:	Monoclonal
Clone ID:	H57-597
Format:	IgG
F/P Ratio:	2-8

Target Details

Gene Name:	Tcrb
Reactivity:	Mouse
Immunogen:	Anti-TCRbeta Antibody (Monoclonal) was produced by repeated immunizations with TCRbeta antigen.
Purity/Specificity:	Fluorescein conjugated TCRbeta Monoclonal Antibody was purified from tissue culture supernatant via affinity chromatography and is directed against mouse TCRbeta. Cross reactivity with TCRbeta from other sources has not been tested. Anti-TCRbeta is conjugated with FITC under optimal conditions and the solution is free of unconjugated FITC.
Relevant Links:	• GeneID - 21577

Application Details

Tested Applications:	FC
Application Note:	Anti-TCR beta is tested for FLOW and useful for Immunofluorescence and Immunoprecipitation using mouse spleen cells, or an appropriate cell type (where indicated). Researchers should determine optimal titers for applications that are not stated.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
FC:	10 μ L/10 ⁶ cells (0.1 μ g)
IF:	User Optimized
IP:	User Optimized

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	0.5mg/mL by UV absorbance at 280 nm
Buffer:	0.01 M Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.09% (w/v) Sodium Azide
Stabilizer:	0.1% Gelatin


Shipping & Handling

Shipping Condition: Wet Ice

Storage Condition: Store vial at 4° C prior to opening. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. DO NOT FREEZE. This product is light sensitive.

Expiration: Expiration date is six (6) months from date of receipt.

Images



Flow Cytometry

Flow Cytometry of anti-TCRbeta Fluorescein Conjugated Monoclonal Antibody. Cells: C57BL/6 splenocytes. Stimulation: none. Antibody: (Black Line) FITC Armenian Hamster IgG isotype control; (Blue Line) Fluorescein Anti-TCRbeta antibody.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.