

Datasheet for 200-B26-N91**CD154 (CD40L) Allophycocyanin Antibody****Overview**

Description:	Anti-CD154 (CD40L) (Armenian Hamster) Allophycocyanin Conjugated Monoclonal Antibody - 200-B26-N91
Item No.:	200-B26-N91
Size:	100 µg
Applications:	FC
Reactivity:	Mouse
Host Species:	Hamster Armenian

Product Details

Background:	CD154 is a 39 kD TNF superfamily member also known as CD40 ligand, gp39, T-BAM, TRAP, and Ly-62. CD154 is an accessory molecule expressed predominantly on activated CD4+ lymphocytes that bind CD40. CD154 plays an important role in T-B cell costimulation. The MR1 antibody has been reported to inhibit the activation of T and B lymphocytes in vitro and antigen-specific lymphocyte responses in vivo.
Synonyms:	CD40 ligand, CD40-L, T-cell antigen Gp39, TNF-related activation protein, TRAP, Tumor necrosis factor ligand superfamily member 5, CD154, Cd40lg, Cd40l, Tnfsf5
Host Species:	Hamster Armenian
Conjugate:	Allophycocyanine (APC)
Clonality:	Monoclonal
Clone ID:	MR1
Format:	IgG
F/P Ratio:	1-2

Target Details

Gene Name:	Cd40lg
Reactivity:	Mouse
Immunogen Type:	Other

Immunogen:	Anti-CD154 Antibody (Monoclonal) was produced by repeated immunizations with Activated mouse Th1 clone D1.6.
Purity/Specificity:	Allophycocyanin conjugated CD154 Monoclonal Antibody was purified from tissue culture supernatant via affinity chromatography and is directed against mouse CD154. Cross reactivity with CD154 from other sources has not been tested. Anti-CD154 is conjugated with APC under optimal conditions and the solution is free of unconjugated APC.
Relevant Links:	• UniProtKB - P27548 • NCBI - NP_035746.2 • GeneID - 21947

Application Details

Tested Applications:	FC
Application Note:	Anti-CD154 is tested for FLOW and useful for Immunohistochemistry using mouse spleen cells, or an appropriate cell type (where indicated). Researchers should determine optimal titers for applications that are not stated.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
FC:	10 μ L/10 ⁶ cells (0.1 μ g)
IHC:	User Optimized

Formulation


Physical State:	Liquid (sterile filtered)
Concentration:	0.2mg/mL Sufficient to run approximately 100 tests
Buffer:	0.01 M Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.09% (w/v) Sodium Azide
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Wet Ice
Storage Condition:	Store vial at 4° C prior to opening. Dilute only prior to immediate use. This product is stable at 4° C as an undiluted liquid. Use subdued lighting during handling and incubation of cells prior to analysis. Store reagent in the dark. DO NOT FREEZE.

Expiration: Expiration date is six (6) months from date of receipt.

Images



Flow Cytometry

Flow Cytometry of anti-CD154 Allophycocyanin Conjugated Monoclonal Antibody. Cells: C57BL/6 mouse splenic T cells. Stimulation: stimulated with PMA + ionomycin for 6 hours. Antibody: (Dotted Line) APC Armenian Hamster IgG isotype control; (Solid Line) Allophycocyanin Anti-CD154 antibody.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.