

Datasheet for 203-401-C41S**Interferon Gamma Antibody****Overview**

Description:	Anti-Chicken Interferon gamma (RABBIT) Antibody - 203-401-C41S
Item No.:	203-401-C41S
Size:	25 µL
Applications:	ELISA, WB
Reactivity:	Chicken
Host Species:	Rabbit

Product Details

Background:	Interferon-gamma (IFN-gamma) is a dimerized soluble cytokine that is the only member of the type II class interferon. This interferon was originally called macrophage-activating factor, a term now used to describe a larger family of proteins to which IFN-gamma belongs. IFN-gamma, or type II interferon, is a cytokine that is critical for innate and adaptive immunity against viral and intracellular bacterial infections and for tumor control. Aberrant IFN-gamma expression is associated with a number of autoinflammatory and autoimmune diseases. The importance of IFN-gamma in the immune system stems in part from its ability to inhibit viral replication directly, but, most important, derives from its immunostimulatory and immunomodulatory effects. IFN-gamma is produced predominantly by natural killer (NK) and natural killer T (NKT) cells as part of the innate immune response, and by CD4 and CD8 cytotoxic T lymphocyte (CTL) effector T cells once antigen-specific immunity develops.
Synonyms:	rabbit anti-Interferon gamma, IFNG, rabbit anti-IFN gamma
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	IFN gamma
Reactivity:	Chicken
Immunogen Type:	Recombinant Protein

Immunogen:	This protein A purified antibody was prepared from whole rabbit serum produced by repeated immunizations with a recombinant protein raised in yeast corresponding to 145 amino acids of mature chicken IFN gamma.
Purity/Specificity:	This product was protein A purified from monospecific antiserum by chromatography. This antibody is specific for chicken IFN gamma protein. A BLAST analysis was used to suggest cross-reactivity with IFN gamma from chicken based on 100% homology; ring-necked pheasant, turkey, quail, and junglefowl based on 91-97% homology with the immunizing sequence. Cross-reactivity with IFN gamma from other sources has not been determined.
Relevant Links:	• UniProtKB - D3GGX1 • NCBI - NP_990480.1 • GenelD - 396054

Application Details

Tested Applications:	ELISA, WB
Application Note:	This protein A purified IFN-gamma antibody has been tested by ELISA and Western blotting. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 16.7 kDa in size corresponding to chicken IFN gamma by western blotting in the appropriate cell lysate or extract.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:20,000
WB:	1:500-1:2,000

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1.24 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Sodium Azide
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Dry Ice
----------------------------	---------

Storage Condition: Store vial at -20° C or below prior to opening. This vial contains a relatively low volume of reagent (25 µL). To minimize loss of volume dilute 1:10 by adding 225 µL of the buffer stated above directly to the vial. Recap, mix thoroughly and briefly centrifuge to collect the volume at the bottom of the vial. Use this intermediate dilution when calculating final dilutions as recommended below. Store the vial at -20°C or below after dilution. Avoid cycles of freezing and thawing.

Expiration: Expiration date is three (3) months from date of receipt.

Images



Western Blot

Western blot using Rockland's protein-A purified anti-chicken IFN gamma antibody shows detection of recombinant chicken IFN gamma at 16.7 kDa, raised in yeast.

Primary antibody was diluted to 1µg/mL. 3% BSA from Rockland's BSA-30 (Bovine Serum Albumin Solution) was used for blocking. Secondary antibody 611-131-122 (Goat anti-Rabbit IgG IRDye 800) was used at 1:20,000.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.